



The Stethoscope Company
... and So Much More

2016 Catalog





Over 40 years ago, Prestige Medical pioneered the concept of pairing diagnostic medical instruments and nursing accessories with medical apparel in retail uniform stores across the United States.

Today, Prestige Medical has built a reputation of excellence with nurses, students, and other healthcare professionals in the U.S. and globally as we have expanded our innovative personalization of healthcare instruments concept to six continents. A relentless focus on stethoscopes, aneroid sphygmomanometers, pulse oximeters, digital blood pressure monitors, medical apparel, medical tools, and nursing accessories makes Prestige Medical the industry leader in outfitting medical professionals.



Cathi Eicher
Director, Uniform Sales



Seth Sternberger
Director, DME Sales



Fay Sicari
Director, College Sales



Joshua Caudillo
Sales



Brian Spangenberg
Sales



Maribel Solano
Manager, Customer Service



Vanessa Valdes
New Accounts Specialist



Andrea Parocua
Receptionist



Carol Bilderbach
Customer Service



Yasmin Aquino
Customer Service



Mercedes Williams
Customer Service



Angela Kiamco
Customer Service



Diagnostics & Instruments

- 2** Stethoscopes
- 11** Veterinary Products
- 24** 3M Littmann®
- 28** Aneroid Sphygmomanometers
- 36** Combination Kits
- 38** Prestige Basics
- 42** Home Healthcare
- 46** Diagnostic Instruments
- 52** EMS Products
- 54** Scissors & Instruments
- 60** Apparel & Accessories



STETHOSCOPES

Why are there so many stethoscopes to choose from and what makes some stethoscopes better than others? Generally speaking, there are three main design features of a stethoscope which impacts its performance.

e Eartips

The eartip is critical to the optimal performance of a stethoscope. It's designed to conform to the outer ear canal to seal out extraneous sound. However, each individual has differences in the size and shape of their ear canals. By following these design requirements, Prestige Medical has developed two sizes of soft silicone eartips for Clinical Series stethoscopes. Silicone provides superior softness, conformability, and resists exposure to chemicals and oils. Prestige Medical also provides soft PVC eartips for most other stethoscopes to enhance their performance and comfort.


t Tubing

Stethoscope tubing is designed to transfer acoustic sound waves channeled by the stethoscope's diaphragm and bell to the practitioner during auscultation. The thickness of the tube is an important factor in enhancing the efficiency of sound transference. Thicker stethoscope tubing transfers sound waves more efficiently while minimizing acoustic interference from the ambient environment. Most Prestige Medical stethoscopes incorporate thick walled tubing to help deliver optimal performance.

m Metal

There are three metals from which most Prestige Medical stethoscopes are constructed. Stainless steel is an alloy prized for its durability and resistance to corrosion, rust, and staining. It is also an excellent metal for the transference of sound waves. Aluminum is very light weight and economical which lends a significant functionality benefit to those who prefer the convenience of a lighter weight instrument. Zinc is excellent for manufacturing complex stethoscope parts economically such as the Sprague-Rappaport chestpiece while still providing very good acoustic characteristics.

Anatomy of a Stethoscope



THE CLINICAL SERIES

Two Piece Box

- ID Tag Included
- Extra Eartips Included
- User Manual & Warranty

The Clinical Series Warranty

All Clinical Series stethoscopes feature a Lifetime Warranty with Free Lifetime Replacement Parts.





Clinical Cardiology®

Model : **128** (Box)

- Enhanced traditional stainless steel cardiology stethoscope
- Includes pediatric diaphragm assembly and interchangeable adult bell
- 2-in-1 tube configuration for optimal sound transference
- Fiberglass diaphragm is used to hear higher frequency heart and respiratory sounds
- Stainless steel bell is used to hear lower frequency heart sounds
- Lifetime Warranty with Free Lifetime Replacement Parts

eartips Soft Silicone

tubing 2-in-1 Thick Walled PVC

metal Stainless Steel

Accessories Included	Two sizes of soft eartips, Replacement diaphragm, Bell and diaphragm assembly
Overall Length	Appx. 27"*
Weight	Appx. 7.0 oz*
Ideal For	Cardiologists Respiratory Specialists General Practitioners Internists ER Medical Students EMTs

* Due to the handmade craftsmanship and material attributes, actual specifications may vary

Tube Colors - Solid:





Clinical Lite™

Models : **121** (Two Piece Box), **S121** (Clamshell)

- A classic dual head stethoscope enhanced by light weight anodized aluminum construction
- Fiberglass diaphragm side is used to hear S1 and S2 heart sounds and Korotkoff sounds
- Bell side is used to hear S3 and S4 heart sounds
- Weighs just 3.9 oz*
- Lifetime Warranty with Free Lifetime Replacement Parts

eartips Soft Silicone

tubing Thick Walled PVC

metal Aluminum

Accessories Included	Two sizes of soft eartips, Replacement diaphragm
Overall Length	Appx. 31" *
Weight	Appx. 3.9 oz*
Ideal For	General Practitioners Nurses Nursing Students Medical Students EMTs

* Due to the handmade craftsmanship and material attributes, actual specifications may vary

Tube Colors



BLK
Black

AQS
Aqua
Sea

HPK
Hot Pink

HUN
Hunter



ICO
Island
Coral

NAV
Navy

N-BLU
Neon
Blue

N-GRN
Neon
Green



N-PNK
Neon
Pink

PAS
Passion

PUR
Purple

RED
Red



ROY
Royal

STE
Stealth

WOR
Wild
Orchid





Clinical I[®]

Models : **126** (Two Piece Box), **S126** (Clamshell)

- A classic dual head stethoscope enhanced with all-stainless steel construction
- Fiberglass diaphragm side is used to hear S1 and S2 heart sounds and Korotkoff sounds
- Bell side is used to hear S3 and S4 heart sounds
- Excellent for all higher level auscultation diagnostics
- Lifetime Warranty with Free Lifetime Replacement Parts

eartips Soft Silicone

tubing Thick Walled PVC


metal Stainless Steel

Accessories Included	Two sizes of soft eartips, Replacement diaphragm
Overall Length	Appx. 31" *
Weight	Appx. 5.8 oz*
Ideal For	General Practitioners Critical Care Nurses Internists Medical Students ER EMTs

** Due to the handmade craftsmanship and material attributes, actual specifications may vary*

Tube Colors

			
BLK Black	AQS Aqua Sea	HPK Hot Pink	HUN Hunter
			
ICO Island Coral	NAV Navy	N-BLU Neon Blue	N-GRN Neon Green
			
N-PNK Neon Pink	PAS Passion	PUR Purple	RED Red
			
ROY Royal	STE Stealth	WOR Wild Orchid	





Clinical I® - Pediatric Edition

Model : **126-PED** (Two Piece Box)

- The Clinical I® stethoscope chestpiece is sized to be acoustically precise for use with children
- Lifetime Warranty with Free Lifetime Replacement Parts

eartips Soft Silicone

tubing Thick Walled PVC

metal Stainless Steel

Accessories Included: Two sizes of soft eartips, replacement diaphragm

Overall Length: Appx. 31"*

Weight: Appx. 4.5 oz*


Ideal For: Pediatricians, General Practitioners, Critical Care Nurses, Internists, Medical Students, ER, EMTs

* Due to the handmade craftsmanship and material attributes, actual specifications may vary

Tube Color



BLK
Black



Clinical I® - Infant Edition

Model : **126-INF** (Two Piece Box)

- The Clinical I® stethoscope chestpiece is sized to be acoustically precise for use with infants
- Lifetime Warranty with Free Lifetime Replacement Parts

eartips Soft Silicone

tubing Thick Walled PVC

metal Stainless Steel

Accessories Included: Two sizes of soft eartips, replacement diaphragm

Overall Length: Appx. 31"*

Weight: Appx. 4.5 oz*


Ideal For: Pediatricians, General Practitioners, Critical Care Nurses, Internists, Medical Students, ER, EMTs

* Due to the handmade craftsmanship and material attributes, actual specifications may vary

Tube Color



BLK
Black





Clinical I® - Gold Edition

Model : **126-G** (Two Piece Box)

- The all-stainless steel construction is enhanced with 22K gold plating
- Lifetime Warranty with Free Lifetime Replacement Parts



Soft Silicone



Thick Walled PVC



22K Gold Plated Stainless Steel

Accessories Included: Two sizes of soft eartips, replacement diaphragm

Overall Length: Appx. 31"*

Weight: Appx. 5.8 oz*


Ideal For: General Practitioners, Critical Care Nurses, Internists,
Medical Students, ER, EMTs

** Due to the handmade craftsmanship and material attributes, actual specifications may vary*

Tube Color



BLK
Black



Clinical I® - Teaching Edition

Model : **126-T** (Two Piece Box)

- The Clinical I® stethoscope is integrated into two binaural headsets to permit dual auscultation between instructor and student
- Lifetime Warranty with Free Lifetime Replacement Parts



Soft Silicone



Thick Walled PVC



Stainless Steel

Accessories Included: Two sizes of soft eartips, replacement diaphragm

Overall Length: Appx. 45"*

Weight: Appx. 8.7 oz*


Ideal For: Instructors, Students

** Due to the handmade craftsmanship and material attributes, actual specifications may vary*

Tube Color



BLK
Black





Clinical Stereo™

Model : 131 (Two Piece Box)

- Revolutionary dual chamber chestpiece and tube design enables true stereo auscultation
- Fiberglass diaphragm is used to hear higher frequency heart and respiratory sounds
- U.S. Patent #D625804
- Lifetime Warranty with Free Lifetime Replacement Parts

e Soft Silicone

t 2-in-1 Thick Walled PVC

m Zinc

Accessories Included: Two sizes of soft eartips, replacement diaphragm

Overall Length: Appx. 27"*

Weight: Appx. 5.4 oz*

Ideal For: Cardiologists, Respiratory Specialists, General Practitioners, Internists, ER

* Due to the handmade craftsmanship and material attributes, actual specifications may vary

Tube Colors



BLK
Black



HPK
Hot Pink




NAV
Navy



PUR
Purple



ROY
Royal



PUR
Purple



Clinical Sprague™

Model : 120 (Two Piece Box)

- Emulates the design of the classic HP Rappaport-Sprague stethoscope
- Specially designed air-filled latex eartips provide extraordinary comfort and sound sealing ability
- Interchangeable bells with form fitted non-chill sleeves for enhanced patient comfort
- Lifetime Warranty with Free Lifetime Replacement Parts

e Air Filled Latex

t Dual Thick Walled PVC

m Zinc

Accessories Included: Firm PVC eartips, adult and pediatric diaphragms, infant and pediatric bells, pediatric diaphragm assembly, disc key

Overall Length: Appx. 24"*

Weight: Appx. 7.3 oz*


Ideal For: General Practitioners, Internists, Nurses, ER, Medical Students, EMTs

* Due to the handmade craftsmanship and material attributes, actual specifications may vary

Tube Color



BLK
Black





Veterinary Clinical I®

Model : **V126** (Box)

- A classic dual head stethoscope enhanced with all stainless steel construction
- Features 2" longer tubing to enable easier auscultation of larger animals
- Lifetime Warranty with Free Lifetime Replacement Parts

e Soft Silicone

t Thick Walled PVC

m Stainless Steel

Accessories Included: Two sizes of soft eartips, replacement diaphragm

Overall Length: Appx. 33"*

Weight: Appx. 5.9 oz*

Ideal For: Veterinarians

* Due to the handmade craftsmanship and material attributes, actual specifications may vary

Tube Colors




BLK
Black



NAV
Navy



PUR
Purple



Veterinary Otoscope

Model : **V3420-L** (Slide Pack)

- Features a long specula
- Twist activated
- Dual-focused standard illumination
- Includes five sizes of reusable specula (cold autoclavable)
- Otoscope head can be removed to convert into penlight
- Includes two AAA batteries



Veterinary Cloisonné Emblem Pins

Each emblem pin is handmade using old world craftsmanship. Painstaking attention to detail and fine gold plating enhance each pin and provide a distinguished appearance. One pin per package.



1049
Veterinary
Medicine



1050
Veterinary
Technician



1051
Registered
Veterinary
Technician

Sprague-Rappaport Stethoscope

Models : **122** (Two Piece Box), **S122** (Clamshell)

- The most widely used stethoscope in the healthcare industry
- 5-in-1 design enables use with adult, pediatric, and infant patients
- Full accessory pouch included
- Lifetime Warranty with Free Lifetime Replacement Parts

e Soft PVC

f Dual Thick Walled PVC

m Zinc

Accessories Included: Soft PVC eartips, adult and pediatric diaphragms, pediatric disc and diaphragm assembly, pediatric and infant bells

Overall Length: Appx. 30"*

Weight: Appx. 7.2 oz*

Ideal For: General Practitioners, Nurses, Nursing / Medical Students, ER, EMTs

** Due to the handmade craftsmanship and material attributes, actual specifications may vary*

Tube Colors - Solids

			
BLK Black	AQS Aqua Sea	BUR Burgundy	CBL Ciel
			
GRY Grey	HPK Hot Pink	HUN Hunter	ICO Island Coral
			
NAV Navy	N-BLU Neon Blue	N-GRN Neon Green	N-PNK Neon Pink
			
PAS Passion	PUR Purple	RED Red	ROY Royal
			
TEA Teal	STE Stealth	WOR Wild Orchid	

Frosted


F-PUR
Purple



Sprague-Rappaport - Gold Edition

Model : **122-G** (Two Piece Box)

- The Sprague-Rappaport stethoscope is enhanced with 22K gold plating
- Full accessory pouch included
- Lifetime Warranty with Free Lifetime Replacement Parts



Soft PVC



Dual Thick Walled PVC



22K Gold Plated Zinc

Accessories Included: Soft PVC eartips, adult and pediatric diaphragms, pediatric disc and diaphragm assembly, pediatric and infant bells

Overall Length: Appx. 30"*

Weight: Appx. 7.2 oz*


Ideal For: General Practitioners, Nurses, Nursing / Medical Students, ER, EMTs

** Due to the handmade craftsmanship and material attributes, actual specifications may vary*

Tube Color



BLK
Black



Sprague-Rappaport - Teaching Edition

Model : **S122-T** (Clamshell)

- Incorporates a standard Sprague-Rappaport chestpiece and two separate tube and binaural sets
- Excellent for instructors to teach students the finer points of auscultation
- Full accessory pouch included
- Lifetime Warranty with Free Lifetime Replacement Parts



Soft PVC



Dual Thick Walled PVC



Zinc

Accessories Included: Soft PVC eartips, adult and pediatric diaphragms, pediatric disc and diaphragm assembly, pediatric and infant bells

Overall Length: Appx. 31.5"*

Weight: Appx. 8.3 oz*

Ideal For: Nursing Instructors, Students

** Due to the handmade craftsmanship and material attributes, actual specifications may vary*

Tube Color



BLK
Black



SpragueLite®

Models : **124** (Two Piece Box), **S124** (Clamshell)

- All the best features of the traditional Sprague-Rappaport, but with significantly less weight
- Versatile 5-in-1 design enables use with adult, pediatric, and infant patients
- Full accessory pouch included
- U.S. Patent #D390660
- Lifetime Warranty with Free Lifetime Replacement Parts

e Soft PVC

t Thick Walled PVC

m Zinc

Accessories Included: Soft PVC eartips, adult and pediatric diaphragms, pediatric disc and diaphragm assembly, pediatric and infant bells

Overall Length: Appx. 32"*

Weight: Appx. 4.9 oz*

Ideal For: General Practitioners, Nurses, Nursing / Medical Students, ER, EMTs

** Due to the handmade craftsmanship and material attributes, actual specifications may vary*

Tube Colors



BLK
Black


HPK
Hot Pink

NAV
Navy



PUR
Purple

ROY
Royal



HPK
Hot Pink

REPLACEMENT PARTS

Prestige Replacement Parts



To view all of our stethoscope accessories and replacement parts, please visit:

www.prestigemedical.com

Dual Head Stethoscope

Models : **108** (Two Piece Box), **S108** (Clamshell)

- Traditional aluminum dual head stethoscope
- Lifetime Warranty with Free Lifetime Replacement Parts

e Soft Silicone **t** Standard PVC **m** Aluminum

Accessories Included: Soft PVC eartips
Overall Length: Appx. 32"*
Weight: Appx. 3.8 oz*
Ideal For: General Practitioners, Nurses, Nursing Students

* Due to the handmade craftsmanship and material attributes, actual specifications may vary

Tube Colors



Ultra-Sensitive Dual Head Stethoscope

Model : **S125** (Clamshell)

- Traditional dual head design incorporating an ultra-sensitive random fiber diaphragm
- Enhanced zinc chestpiece and thick walled tubing
- Lifetime Warranty with Free Lifetime Replacement Parts

e Soft Silicone **t** Thick Walled PVC **m** Zinc

Accessories Included: Soft PVC eartips
Overall Length: Appx. 32"*
Weight: Appx. 5.1 oz*
Ideal For: General Practitioners, Nurses, Nursing / Medical Students

* Due to the handmade craftsmanship and material attributes, actual specifications may vary

Tube Colors



Dual Head - Pediatric Edition

Model : **S108-P** (Clamshell)

- The traditional dual head stethoscope is sized to be acoustically precise for diagnostic use with children
- Lifetime Warranty with Free Lifetime Replacement Parts

 Soft Silicone

 Standard PVC

 Aluminum

Accessories Included: Soft PVC eartips

Overall Length: Appx. 31.5"*

Weight: Appx. 3.5 oz*

Ideal For: Pediatricians

** Due to the handmade craftsmanship and material attributes, actual specifications may vary*

Tube Colors



BLK
Black

HPK
Hot Pink



HPK
Hot Pink

Dual Head - Infant Edition

Model : **S108-I** (Clamshell)

- The traditional dual head stethoscope is sized to be acoustically precise for diagnostic use with Infants
- Lifetime Warranty with Free Lifetime Replacement Parts

 Soft Silicone

 Standard PVC

 Aluminum

Accessories Included: Soft PVC eartips

Overall Length: Appx. 31.5"*

Weight: Appx. 3.5 oz*

Ideal For: Pediatricians


** Due to the handmade craftsmanship and material attributes, actual specifications may vary*

Tube Colors



BLK
Black

HPK
Hot Pink



BLK
Black

Clear Sound™

Model : **S107** (Clamshell)

- Ultra-sensitive acoustics
- Accurate assessment through up to 10 layers of folded towels
- Chestpiece made from impact-resistant acrylic resin
- Lifetime Warranty with Free Lifetime Replacement Parts



Soft PVC



Thick Walled PVC



Acrylic Resin

Accessories Included: Soft PVC eartips, ID tag

Overall Length: Appx. 30"*

Weight: Appx. 4.2 oz*

Ideal For: Nurses

** Due to the handmade craftsmanship and material attributes, actual specifications may vary*

Tube Colors - Solid



BLK
Black



HPK
Hot Pink




PLU
Plum

Frosted



F-PUR
Purple



PLU
Plum

Clear Sound™ - Breast Cancer Awareness Edition

Model : **S107-PR*** (Clamshell)

- All the same features as the standard Clear Sound™
- Unique Pink Ribbon imprinted inside the chestpiece
- A portion of the sales of this stethoscope is donated to the American Cancer Society
- Lifetime Warranty with Free Lifetime Replacement Parts



Soft PVC



Thick Walled PVC



Acrylic Resin

Accessories Included: Soft PVC eartips, ID tag

Overall Length: Appx. 30"*

Weight: Appx. 4.2 oz*

Ideal For: Nurses

** Due to the handmade craftsmanship and material attributes, actual specifications may vary*



Clear Sound™ - Heart Edition

Model : **S107-H** (Clamshell)

- All the same features as the standard Clear Sound™
- Unique heart-shaped chestpiece
- Patent #D566837
- Lifetime Warranty with Free Lifetime Replacement Parts

 Soft PVC

 Thick Walled PVC

 Acrylic Resin

Accessories Included: Soft PVC eartips, ID tag

Overall Length: Appx. 31"*

Weight: Appx. 4.5 oz*

Ideal For: Nurses

** Due to the handmade craftsmanship and material attributes, actual specifications may vary*

Tube Color - Solid



Red with Black

Tube Colors - Frosted



F-PUR
Purple



Clear Sound™ - Printed Edition

Model : **S107** (Clamshell)

- All the same features as the standard Clear Sound™
- Designs imprinted inside the chestpiece
- Lifetime Warranty with Free Lifetime Replacement Parts

 Soft PVC

 Thick Walled PVC

 Acrylic Resin

Accessories Included: Soft PVC eartips, ID tag

Overall Length: Appx. 30"*

Weight: Appx. 4.2 oz*

Ideal For: Nurses

** Due to the handmade craftsmanship and material attributes, actual specifications may vary*


Black Tubes - Prints



SML
Smiley Face



LUV
Betty Boop®
Luv-A-Nurse



Single Head Stethoscope

Models : **106** (Two Piece Box), **S106** (Clamshell)

- Ergonomically designed lightweight aluminum chestpiece
- Designed to fit under the cuff during blood pressure assessment or for basic diagnostics
- Lifetime Warranty with Free Lifetime Replacement Parts

e Soft Silicone **t** Standard PVC **m** Aluminum

Accessories Included: Soft PVC eartips

Overall Length: Appx. 31.5"*

Weight: Appx. 2.8 oz*

Ideal For: Nurses, Nursing Students, Blood Pressure Measurement

** Due to the handmade craftsmanship and material attributes, actual specifications may vary*

Tube Colors



BLK
Black



Single Patient Stethoscope

Model : **100** (Poly Bag)

- Disposable stethoscope designed for single patient use
- Prevents cross-contamination
- Yellow tubing only

e Firm PVC **t** Standard PVC **m** Plastic

Accessories Included: None

Overall Length: Appx. 30"*

Weight: Appx. 1.9 oz*

Ideal For: General Assessment in cross-contamination hazard or quarantine conditions

** Due to the handmade craftsmanship and material attributes, actual specifications may vary*



Small Silicone Clinical Eartips

Soft premium push-on silicone eartips. Small size.

Black only.



CE-SIL-S - Unpackaged
S-CE-SIL-S - Slide Pack

Fits:

S121, 121
S126, 126
127
128

Threaded Eartips

Threaded soft PVC eartips. Two sizes.

Fits:

112
S122, 122
S124, 124



BLK
Black

GRY
Grey

DET-122-S - Small size - Unpackaged

S-DET-122-S - Small size - Slide pack

DET-122-L - Large size - Unpackaged

S-DET-122-L - Large size - Slide pack

Large Silicone Clinical Eartips

Soft premium push-on silicone eartips.

Large size.

Black only.



CE-SIL-L - Unpackaged
S-CE-SIL-L - Slide Pack

Fits:

S121, 121
S126, 126
127
128

Clear PVC Clinical Eartips

Soft premium push-on clear PVC eartips.



Fits:

S121, 121
S126, 126
127
128

CE-PVC - Unpackaged
S-CE-PVC - Slide Pack

Mushroom Eartips

Model : **S-SET**

Soft PVC eartips.
Slide pack.

Black only.



Fits:

103
104
S106, 106
S108, 108
S125

Push-On Eartips

Soft push-on PVC eartips. Two sizes.



BLK
Black

GRY
Grey

DET-SN-S - Small size - Unpackaged
S-DET-SN-S - Small size - Slide pack
DET-SN-L - Large size - Unpackaged
S-DET-SN-L - Large size - Slide pack

Fits:

S107
120
S121, 121
S126, 126
127
128
131

Soft Premium Eartips

Air-filled latex eartips.

Black only.



Fits:

120
S122, 122
S124, 124
131

DET-120 - Unpackaged
DET-S120 - Slide Pack

Firm Plastic Eartips

Model : **ET**

Firm PVC eartips.
Slide pack.

White only.



Fits:

103
104

Stethoscope Tape Holder

Model : **754** (Slide Pack)


- Attaches to any stethoscope
- Keeps a roll of 1" medical tape handy
- A perfect counter-weight to help balance the stethoscope around the user's neck
- Made in the U.S.A.



Quick Equip™ Stethoscope Light

Model : **S205** (Slide Pack)

An LED light that snaps onto a stethoscope chestpiece. Compatible with Prestige models 121/S121, 126/S126, 124/S124, 108/S108, S125. *Not recommended for eye examinations.*



Stethoscope Watch

Model : **1688** (Slide Pack)


Back clip attaches to most stethoscope tubing. 24 hour face with step second hands. Quartz movement and a water resistant case. Includes an LR626 battery.



Medical Symbols Stethoscope Watch

Model : **1689** (Slide Pack)


Back clip attaches to most stethoscope tubing. 24 hour face with step second hands. Quartz movement and a water resistant case. Includes an LR626 battery.



Two-Sided ID Tag

Model : **S4** (Slide Pack)

Double-sided ID tag enables identification to be written on both sides. Universal fit.



Sprague ID Tag

Model : **S7-CLR** (Slide Pack)

Single-sided ID Tag. Designed for the Sprague-Rappaport or dual tube stethoscope.



Printed Stethoscope ID Tag

Model : **S8** (Slide Pack)


Single-sided stethoscope ID tags featuring a "domed" printed face. A virtually indestructible adhesive closure enables secure attachment to any stethoscope tube.



Hip Clip™ Stethoscope Holder

Model : **755** (Poly Bag)

A convenient stethoscope holder that affixes to the user's belt or waist band. Eliminates the fatigue of wearing a stethoscope around the neck.




CharMED™ Crystal Stethoscope Charms

Model : **CEG-0001** (Header Card)


- The perfect way to create a sense of uniqueness and style for the health-care professional
- Fully adjustable and will fit most stethoscope or aneroid sphygmomanometer tubes
- Made of an antibacterial zinc alloy



AQM
Aquamarine



BYL
Bright Yellow



CHY
Cherry



FUS
Fuschia



ICE
Ice



IRD
Iridescent



LGS
Lemongrass



OMB
Ombre



ROY
Royal



VIO
Violet



PAW
Paw Print



SML
Smiley Face



CNA
Certified Nursing
Assistant



LPN
Licensed Practical
Nurse



MAS
Medical
Assistant



RNB
Registered
Nurse



PAT
Physicians
Assistant



PKB
Pink Ribbon*



HRT
Heart



BFT
Baby Feet



MTC
Maltese Cross



SOL
Star of Life



EKB
EKG on Black



LDB
Lady Bug


REPLACEMENT PARTS

Prestige Replacement Parts



To view all of our stethoscope accessories and replacement parts, please visit:
www.prestigemedical.com

3200 ELECTRONIC



Burgundy

The Electronic Model 3200 goes even further than the Model 3100 - with on-board recording and playback capabilities. Also, Bluetooth® technology lets you transmit sounds real-time to your PC, which can then be further analyzed, attached to medical records, or reviewed online with colleagues. Free heart & lung sound visualization software, and available software that identifies suspected heart murmurs with a Class I indication. It all adds up to an exceptional listening experience compared to acoustic scopes.


Black - 3200BK27

Burgundy - 3200BU

Navy - 3200NB



3100 ELECTRONIC



Navy

This model's exceptional sound quality reduces background noise an average of 85% with its Ambient Noise Reduction Technology and amplifies sounds up to 24 times to hear faint sounds you might otherwise miss. You can auscultate immediately, thanks to the instant-on technology. The power button quickly activates backlight for low-light conditions and the stethoscope stays in battery-saving "sleep mode" until you start your next auscultation. The backlit LCD display includes frequency selection, sound level, remaining battery life and your patient's heart rate.

Black - 3100BK27

Burgundy - 3100BU

Navy - 3100NB

MASTER CARDIOLOGY



Smoke /
Black

Provides the best acoustics Littmann stethoscopes offer. Solid stainless steel single sided chestpiece. Two tubes in one design reduces noise from two tubes rubbing together. Soft sealing eartips. A patented special procedure adapter allows auscultation of difficult body parts, and pediatric or emaciated patients.

Black (22") - 2159

Plum - 2167

Black (27") - 2160

Burgundy - 2163

Hunter - 2165

Caribbean - 2178

Navy - 2164

All Black - 2161

Brass Finish/Black Tube - 2175

Brass Finish/Hunter Tube - 2183

Smoke Finish/Black Tube - 2176

Smoke Finish/Dark Olive Tube - 2182

CARDIOLOGY S.T.C.



Plum

Innovative design and outstanding acoustics are combined with a unique and easy to grasp Soft-Touch Chestpiece that allows easy movement between auscultation sites. Single sided stainless steel chestpiece provides both bell and diaphragm frequency responses. Diaphragm allows user to alternate between low and high frequency by applying pressure.

Black (27") - 4471

Burgundy - 4472

Navy - 4473

Plum - 4475

CARDIOLOGY III




Copper Finish

Superior versatility and outstanding acoustics with two tunable diaphragms on a stainless steel two-sided chestpiece. Two tubes-in-one design.

Black (22") - 3127	Gray - 3136
Black (27") - 3128	Chocolate - 3137
Burgundy - 3129	Caribbean - 3138
Ceil Blue - 3146	Coral Pink - 3149
Hunter - 3134	Red - 3140
Navy - 3130	Raspberry - 3148
Plum - 3135	All Black - 3131BE
Brass Finish/Black Tube - 3128BRS	
Copper Finish/Chocolate Tube - 3137CPR	
Rainbow Finish/Black Tube - 3152RBW	
Rainbow Finish/Lavender Tube - 3158	
Smoke Finish/Black Tube - 3157SM	
Smoke Finish/Dark Olive Tube - 3166	
Smoke Finish/Sand Tube - 3165	
Smoke Finish/Lemon-Lime Tube - 3161	
Smoke Finish/Turquoise Tube - 3159	

MASTER CLASSIC II




Smoke / Sand

Offers efficiency in stethoscope auscultation that is unprecedented in a general examination stethoscope. The bell and diaphragm modes are incorporated into a single-sided Master chestpiece. No need to remove the chestpiece from the patient and mechanically switch from bell to diaphragm sides.

Black - 2144L	Hunter - 2632
Purple - 2143	Navy - 2147
Burgundy - 2146	Pine Green - 2634
Caribbean - 2630	Ceil Blue - 2633
All Black - 2141	
Smoke Finish/Sand Tube - 2636	


MASTER CLASSIC II TEACHING



The 3M Littmann® Master Classic II Stethoscope has two headsets, which allow instructor and student to perform simultaneous auscultation. Features a single-sided chestpiece with the tunable diaphragm.

Black/Burgundy - 2139

MASTER CLASSIC II VETERINARY




Navy

This version of the Master Classic II features longer tubing for Veterinary professionals. Offers efficiency in stethoscope auscultation that is unprecedented in a general examination stethoscope. The bell and diaphragm modes are incorporated into a single-sided Master chestpiece.

Navy - 1392

CLASSIC III



Orange

Littmann's latest addition to their line. Features a two-sided chestpiece, dual tunable diaphragms, convertible open bell, and improved tubing.

Black - 5620	Burgundy - 5627
Caribbean - 5623	Ciel Blue - 5630
Emerald - 5840	Gray - 5621
Hunter - 5624	Lavender - 5832
Lemon Lime - 5839	Lime Green - 5829
Navy - 5622	Orange - 5629
Pearl Pink - 5633	Plum - 5831
Raspberry - 5626	Turquoise - 5835
Black Finish/Black Tube - 5803	
Smoke Finish/Black Tube - 5811	
Rainbow Finish/Caribbean Tube - 5807	
Copper Finish/Chocolate Tube - 5809	
Smoke Finish/Dark Olive Tube - 5812	
Rainbow Finish/Raspberry Tube - 5806	

CLASSIC II SE



Raspberry

The 3M™ Littmann® Classic II S.E. Stethoscope is the standard of the industry and features 3M's patented tunable diaphragm for improved low-frequency response and the ability to hear both low and high-frequency sounds without turning over the chestpiece. Also features a nonchill rim, single-lumen tubing and solid stainless steel combination chestpiece.

Black - 2201
Ciel Blue - 2813
Navy - 2205
Peach - 2822
Pearl Pink - 2817
Raspberry - 2210
All Black - 2218BE

CLASSIC II SE PEDIATRIC



Rainbow / Caribbean

The Classic II SE Pediatric features a dual frequency, tunable diaphragm that allows the user to listen to both low and high frequency. The diaphragm on one side augments the traditional bell function for the opposite side. Firm pressure allows the user to hear high frequency sounds and light pressure lets them hear low frequency sounds. The specially designed bell and diaphragm chestpiece that is designed, sized and acoustically precise for children.

Black - 2113	Raspberry - 2122
Caribbean - 2119	Red - 2113R
Navy - 2123	Pine Green - 2131
Peach - 2154	Orange - 2155
Royal Blue - 2136	Pearl Pink - 2115
Rainbow Finish/Caribbean Tube - 2153	

CLASSIC II SE INFANT



Royal Blue


The Classic II SE Infant features a dual frequency, tunable diaphragm that allows the user to listen to both low and high frequency. The diaphragm on one side augments the traditional bell function for the opposite side. Firm pressure allows the user to hear high frequency sounds and light pressure lets them hear low frequency sounds. The specially designed bell and diaphragm chestpiece that is designed, sized and acoustically precise for infants.

Black - 2114	Raspberry - 2125
Caribbean - 2124	Red - 2114R
Pearl Pink - 2120	Royal Blue - 2156
Peach - 2158	Orange - 2179
Rainbow Finish/Raspberry Tube - 2157	

CLASSIC II SE TEACHING

The 3M Littmann® Classic II S.E. Stethoscope has two headsets, which allow instructor and student to perform simultaneous auscultation. Single-sided chestpiece with the tunable diaphragm.


Black/Burgundy - 2138



LIGHTWEIGHT II SE

The perfect entry level and affordable choice stethoscope. The tunable diaphragm assures acoustic sensitivity and alternates between high and low frequency sounds. Special lightweight construction.

Ceil Blue




Black - 2450
Burgundy - 2451
Caribbean - 2452
Lilac - 2453
Ceil Blue - 2454
Pearl Pink - 2456
Seafoam - 2455

SELECT

A convenient tunable diaphragm in a single-sided, coated alloy chestpiece providing reliable and clear acoustics. The Select Stethoscope features a tunable diaphragm that allows the user to conveniently alternate between bell and diaphragm modes by applying simple pressure on the chestpiece. Its clear acoustics offer an advantage to healthcare professionals when listening to heart, lung and blood pressure sounds.

Ocean Blue



Black - 2290	Royal - 2298
Burgundy - 2293	Purple - 2294
Caribbean - 2291	Pine Green - 2305
Ceil Blue - 2301	Ocean Blue - 2306
Pearl Pink - 2292	Peach - 2310
Raspberry - 2296	

LITTMANN® PARTS & MATERIALS



Please see page 7 of our 2016 Price List for a complete list of Littmann® replacement parts and materials.

LITTMANN® POINT-OF-PURCHASE MATERIALS

Please see page 8 of our 2016 Price List for a complete list of Littmann® point-of-purchase materials.





ANEROID SPHYGMOMANOMETERS

What makes Prestige Medical a leader in aneroid sphygmomanometers and combination kits ?

■ Lifetime Warranty

Accuracy and performance are what Prestige Medical's aneroid sphygmomanometers are renowned for. Each manometer gauge is individually tested to meet or exceed the stringent guidelines set by the American National Standards Institute.

It is because of this rigorous testing that most Prestige Medical aneroid sphygmomanometer models are backed with a Lifetime Calibration Warranty.

■ More Choices

Offering the largest selection of solid colors and fashion prints in the industry, Prestige Medical is committed to personalizing the instruments of healthcare.

New prints and styles are added based on feedback from our dealers and consumers. Visit us online at www.prestigemedical.com to view all of our styles in each of their unique color offerings by simply clicking on the associated color swatches.



■ Custom Combination Kits

Perfect for group orders.

If you don't see a combination kit here that meets your needs, simply create your own. Prestige Medical offers hundreds of possible combinations.

To create your own custom combination kit, please contact customer service or your sales representative.

Anatomy of an Aneroid Sphygmomanometer



*A gauge is still fully calibrated if the needle is inside or on the line of the calibration rectangle.

Premium Aneroid Sphygmomanometer

Models : **82 (Box)**, **S82 (Clamshell)**

- Our most popular adult size aneroid sphygmomanometer
- Lifetime Calibration Warranty

Features: Latex free, nylon cuff, Index and Range markings, artery indicator mark, gauge holder, owner ID label, standard inflation bulb and air release valve

Cuff Range: 10"–16"*

*Range refers to the arm circumference that the cuff is designed to operate on

Cuff Colors



Cuff Prints



Premium Pediatric & Infant Aneroid Sphygmomanometers

Models : **82-PED (Pediatric)**, **82-INF (Infant)**

- Premium aneroid sphygmomanometer sized for use on infants, children, or small adults
- Lifetime Calibration Warranty

Features: Latex free, nylon cuff, Index and Range markings, artery indicator mark, ID label, gauge holder, standard inflation bulb and air release valve


Cuff Range: Pediatric : 7.25"–10.5"*, Infant : 5"–7.25"*

*Range refers to the arm circumference that the cuff is designed to operate on

Cuff Color



Cuff Print



Premium Large Adult Aneroid Sphygmomanometer

Model : **82-OB** (Box)

- Premium aneroid sphygmomanometer sized for use on large adults
- Lifetime Calibration Warranty

Features: Latex free, nylon cuff, Index and Range markings, artery indicator mark, ID label, gauge holder, standard inflation bulb and air release valve

Cuff Range: 13.5" - 20"*

**Range refers to the arm circumference that the cuff is designed to operate on*

Cuff Colors



Cuff Print



CGS
Chocolate
Golden Swirls



Standard Aneroid Sphygmomanometer

Models : **79** (Box), **S79** (Clamshell)


- Standard adult size aneroid sphygmomanometer
- 10 Year Calibration Warranty

Features: Latex free, nylon cuff, Index and Range markings, artery indicator mark, gauge holder, owner ID label, standard inflation bulb and air release valve

Cuff Range: 10"-16"*

**Range refers to the arm circumference that the cuff is designed to operate on*

Cuff Colors



Cuff Print



CGS
Chocolate
Golden Swirls



Premium Aneroid Sphygmomanometer with Carry Case

Model : **882** (Clear Box)

- Our most popular aneroid sphygmomanometer (82/S82) color-matched with an extra-large 6" x 9" nylon carry case (745)
- Lifetime Calibration Warranty


Features: Latex free, nylon cuff, Index and Range markings, artery indicator mark, gauge holder, owner ID label, standard inflation bulb and air release valve

Cuff Range: 10"–16"*

*Range refers to the arm circumference that the cuff is designed to operate on



Cuff Colors



Cuff Prints



REPLACEMENT PARTS

Prestige Replacement Parts



To view our aneroid sphygmomanometer replacement parts, please visit:

www.prestigemedical.com

3-in-1 Aneroid Sphygmomanometer Set with Carry Case

Model : **882-COM** (Clear Box)

- Features three interchangeable cuff sizes - pediatric, adult, and large adult
- Includes an extra-large 6" x 9" color-matched nylon carry case
- Lifetime Calibration Warranty

Features: Latex free, nylon cuff, Index and Range markings, artery indicator mark, gauge holder, owner ID label, standard inflation bulb and air release valve

Cuff Ranges: Adult : 10" - 16"*
Pediatric : 7.25" - 10.5"*
Large Adult : 13.5" - 20"*

*Range refers to the arm circumference that the cuff is designed to operate on

Cuff Color



BLK
Black

Cuff Print



CGS
Chocolate
Golden Swirls



Premium Extra-Large Adult Aneroid Sphygmomanometer with Carry Case

Model : **882-TH** (Clear Box)

- Premium aneroid sphygmomanometer sized for use on extra-large adults or a thigh
- Lifetime Calibration Warranty

Features: Latex free, nylon cuff, Index and Range markings, artery indicator mark, gauge holder, owner ID label, standard inflation bulb and air release valve, extra-large 6" x 9" nylon carry case

Cuff Range: 16" - 26"*

*Range refers to the arm circumference that the cuff is designed to operate on

Cuff Color



BLK
Black

Cuff Print



CGS
Chocolate
Golden Swirls



Easy-ID™ Latex Free Aneroid Sphygmomanometer

- Unique terra-cotta color of the neoprene inflation system enables easy identification of its latex free properties
- Lifetime Calibration Warranty
- Black only

Features: Latex free, nylon cuff, Index and Range markings, artery indicator mark, gauge holder, owner ID label, standard inflation bulb and air release valve

Cuff Range: Adult : 10"–16"*

Pediatric : 7.25" - 10.5"*

Large Adult : 13.5" - 20"*

*Range refers to the arm circumference that the cuff is designed to operate on

S72-BLK	Adult (Clamshell)
S72-PED-BLK	Pediatric (Clamshell)
S72-OB-BLK	Large Adult (Clamshell)



Cotton Cuff Aneroid Sphygmomanometer

- The soft comfort of a cotton cuff is what sets this model apart
- Lifetime Calibration Warranty
- Black only

Features: Latex free, cotton cuff, Index and Range markings, artery indicator mark, gauge holder, owner ID label, deluxe inflation bulb and air release valve

Cuff Range: Adult : 10"–16"*

Pediatric : 7.25" - 10.5"*

Large Adult : 13.5" - 20"*

*Range refers to the arm circumference that the cuff is designed to operate on

80	Adult (Box)
80-PED	Pediatric (Box)
80-OB	Large Adult (Box)



Single-Hand Dual Tube Aneroid Sphygmomanometer

Model : **883** (Clear Box)

- High performance integrated inflation bulb and manometer gauge for single hand operation
- Dual tube manometer gauge design enables it to be interchangeable with most dual tube cuffs
- Includes an extra-large 6" x 9" color-matched nylon carry case
- Lifetime Calibration Warranty

Features: Latex free, nylon cuff, Index and Range markings, artery indicator mark, gauge holder, owner ID label

Cuff Range: 10" - 16"*

**Range refers to the arm circumference that the cuff is designed to operate on*

Cuff and Case Colors



BLK
Black

NAV
Navy



Single-Hand One Tube Aneroid Sphygmomanometer

Model : **887** (Clear Box)

- High performance integrated inflation bulb and manometer gauge for single hand operation
- Interchangeable only with other one tube cuffs
- Includes an extra-large 6" x 9" color-matched nylon carry case
- Lifetime Calibration Warranty

Features: Latex free, nylon cuff, Index and Range markings, artery indicator mark, gauge holder, owner ID label

Cuff Range: 10" - 16"*

**Range refers to the arm circumference that the cuff is designed to operate on*

Cuff and Case Colors



BLK
Black

NAV
Navy



Aneroid Sphygmomanometer / Sprague-Rappaport Kit

Model : **A2** (Clear Box)

- Our most popular aneroid sphygmomanometer matched with the Sprague-Rappaport stethoscope
- Includes an extra-large 6" x 9" color-matched nylon carry case
- Lifetime Calibration Warranty

Features: Latex free, nylon cuff, Index and Range markings, artery indicator mark, gauge holder, owner ID label, standard inflation bulb and air release valve, full accessory pouch for Sprague-Rappaport

Cuff Range: 10"-16"*

**Range refers to the arm circumference that the cuff is designed to operate on*

Build your own custom kits.
See page 87 for more info.



WAB
Woody
Animals Black

Cuff Colors / Stethoscope Colors



BLK
Black

AQS
Aqua Sea

BUR
Burgundy

CBL
Ciel



HUN
Hunter

ICO
Island Coral

NAV
Navy

N-BLU
Neon Blue
(Black Cuff)



N-GRN
Neon Green
(Black Cuff)

N-PNK
Neon Pink
(Black Cuff)

PAS
Passion

PUR
Purple



ROY
Royal

TEA
Teal

WOR
Wild
Orchid

STE
Stealth
(Black Cuff)

Cuff Prints / Stethoscope Colors



CMG
Green
Camouflage




CMP
Pink
Camouflage



FSB
Floral Stitch
Black



FSO
Floral Stitch
Orchid



LPC
Leopard Print
Cream



LPG
Leopard Print
Grey



LVC
Leaves
Cream



LVG
Leaves
Grey



PBB
Pink Ribbons
& Butterflies*



RHB
Ribbons and
Hearts Black*




RPK
Ribbons and
Hearts Pink*



RPU
Ribbons and
Hearts Purple*



WAB
Woody
Animals Black



WAC
Woody
Animals Cream

Clinical Lite™ Combination Kit

Model : **A121** (Clear Box)

- Our most popular aneroid sphygmomanometer matched with a Clinical Lite™ stethoscope
- Includes an extra-large 6" x 9" color-matched nylon carry case
- Lifetime Calibration Warranty

Features: Latex free, nylon cuff, Index and Range markings, artery indicator mark, gauge holder, owner ID label, standard inflation bulb and air release valve, two sizes of soft eartips, extra diaphragm, stethoscope ID tag

Cuff Range: 10"–16"*

**Range refers to the arm circumference that the cuff is designed to operate on*

Build your own custom kits.
See page 87 for more info.



Cuff Colors / Stethoscope Colors



BLK
Black

AQS
Aqua Sea

HUN
Hunter

ICO
Island Coral



NAV
Navy

N-BLU
Neon Blue
(Black Cuff)

N-GRN
Neon Green
(Black Cuff)

N-PNK
Neon Pink
(Black Cuff)



PAS
Passion

PUR
Purple

ROY
Royal

WOR
Wild Orchid



STE
Stealth
(Black Cuff)

Cuff Prints / Stethoscope Colors



CMG
Green Camouflage




CMP
Pink Camouflage



FSB
Floral Stitch Black



FSO
Floral Stitch Orchid



LPC
Leopard Print Cream



LPG
Leopard Print Grey



LVC
Leaves Cream



LVG
Leaves Grey



PBB
Pink Ribbons & Butterflies*



RHB
Ribbons and Hearts Black*




RPK
Ribbons and Hearts Pink*



RPU
Ribbons and Hearts Purple*



WAB
Woody Animals Black



WAC
Woody Animals Cream

Basic Classic Stethoscope

Model : **117-BLK** (Box)

- Economical Classic stethoscope
- Black only

Accessories Included: Two sizes of soft silicone eartips, replacement diaphragm

Overall Length: Appx. 27"

Weight: Appx. 7.0 oz*

Warranty: One Year

** Due to the handmade craftsmanship and material attributes, actual specifications may vary*



Basic Sprague-Rappaport Stethoscope

Model : **105-BLK** (Box)

- Economical Sprague-Rappaport stethoscope
- Black only

Accessories Included: Soft PVC eartips, adult and pediatric diaphragms, pediatric diaphragm assembly, pediatric and infant bells

Overall Length: Appx. 30"

Weight: Appx. 7.2 oz*

Warranty: One Year

** Due to the handmade craftsmanship and material attributes, actual specifications may vary*



Basic SpragueLite® Stethoscope

Model : **110-BLK** (Box)

- Economical SpragueLite® stethoscope
- Black only

Accessories Included: Soft PVC eartips, adult and pediatric diaphragms, pediatric diaphragm assembly, pediatric and infant bells

Overall Length: Appx. 32"

Weight: Appx. 4.9 oz*

Warranty: One Year

** Due to the handmade craftsmanship and material attributes, actual specifications may vary*



Basic Dual Head Stethoscope

Model : **104-BLK** (Box)

- Economical dual head stethoscope
- Black only

Accessories Included: Soft PVC eartips
Overall Length: Appx. 32"*
Weight: Appx. 3.8 oz*
Warranty: One Year

** Due to the handmade craftsmanship and material attributes, actual specifications may vary*



Basic Single Head Stethoscope

Model : **103-BLK** (Box)

- Economical single head stethoscope
- Black only

Overall Length: Appx. 31.5"*
Weight: Appx. 2.8 oz*
Warranty: One Year

** Due to the handmade craftsmanship and material attributes, actual specifications may vary*



Basic Aneroid Sphygmomanometer

- Economical aneroid sphygmomanometer
- Available in both adult and large adult sizes

Features: Latex free, nylon cuff, Index and Range markings, artery indicator mark, gauge holder, owner ID label, standard inflation bulb and air release valve
Cuff Range: Adult : 10"-16"*, Large Adult : 13.5" - 20"*
Warranty: One Year

**Range refers to the arm circumference that the cuff is designed to operate on*

70	Adult
70-OB	Large Adult

Cuff Color



BLK
Black

Cuff Print



CGS
Chocolate
Golden Swirls



Basic Aneroid Sphygmomanometer / Single Head Kit

Model : **A1-103** (Box)

- Economical aneroid sphygmomanometer / single head combination kit

Features: Latex free, nylon cuff, Index and Range markings, artery indicator mark, gauge holder, owner ID label, standard inflation bulb and air release valve, extra-large 6" x 9" color-matched nylon carry case

Cuff Range: 10"-16"*

Warranty: One Year

* Range refers to the arm circumference that the cuff is designed to operate on

**Cuff Color /
Stethoscope
Color**



BLK
Black

**Cuff Print /
Stethoscope
Color**



CGS
Chocolate
Golden Swirls



Basic Aneroid Sphygmomanometer / Dual Head Kit

Model : **A1-104** (Box)

- Economical aneroid sphygmomanometer / dual head combination kit

Features: Latex free, nylon cuff, Index and Range markings, artery indicator mark, gauge holder, owner ID label, standard inflation bulb and air release valve, extra-large 6" x 9" color-matched nylon carry case

Cuff Range: 10"-16"*

Warranty: One Year

* Range refers to the arm circumference that the cuff is designed to operate on

**Cuff Color /
Stethoscope
Color**



BLK
Black

**Cuff Print /
Stethoscope
Color**



CGS
Chocolate
Golden Swirls



Basic Aneroid Sphygmomanometer / Sprague-Rappaport Kit

Model : **A1-105** (Box)

- Economical aneroid sphygmomanometer / Sprague-Rappaport combination kit

Features: Latex free, nylon cuff, Index and Range markings, artery indicator mark, gauge holder, owner ID label, standard inflation bulb and air release valve, extra-large 6" x 9" color-matched nylon carry case, full accessory pouch for Sprague-Rappaport

Cuff Range: 10"-16"

Warranty: One Year

* Range refers to the arm circumference that the cuff is designed to operate on

**Cuff Color /
Stethoscope
Color**



BLK
Black

**Cuff Print /
Stethoscope
Color**



CGS
Chocolate
Golden Swirls



Basic Aneroid Sphygmomanometer / SpragueLite® Kit

Model : **A1-110** (Box)

- Economical aneroid sphygmomanometer / SpragueLite® combination kit

Features: Latex free, nylon cuff, Index and Range markings, artery indicator mark, gauge holder, owner ID label, standard inflation bulb and air release valve, extra-large 6" x 9" color-matched nylon carry case, full accessory pouch for SpragueLite®

Cuff Range: 10"-16"

Warranty: One Year

* Range refers to the arm circumference that the cuff is designed to operate on

**Cuff Color /
Stethoscope
Color**



BLK
Black

**Cuff Print /
Stethoscope
Color**



CGS
Chocolate
Golden Swirls



WristMate™ Digital Blood Pressure Monitor

Model : **HM-50*** (Box)

- Automatic blood pressure and pulse rate monitor that measures from the wrist
- One-touch operation and clinically tested for accuracy
- Recalls up to 40 prior readings for up to three users
- Includes World Health Organization (WHO) indicator for blood pressure classification

Includes : Two AAA batteries and plastic storage case

Cuff Range: 5.3" - 7.5"*

Warranty: Two Year

**Range refers to the wrist circumference that the cuff is designed to operate on*



HealthMate® Digital Blood Pressure Monitor

Model : **HM-30*** (Box)

- Fully automatic blood pressure and pulse rate monitor
- One-touch operation and clinically tested for accuracy
- Recalls up to 99 prior readings
- World Health Organization (WHO) indicator for BP classification

Includes : Four AA batteries, DC power port and nylon storage case

Cuff Range: Adult: 9" - 13"*

Warranty: Two Year

**Range refers to the arm circumference that the cuff is designed to operate on*



HealthMate® Digital Blood Pressure Monitor Large Adult Edition

Model : **HM-30-OB*** (Box)

- Fully automatic blood pressure and pulse rate monitor
- One-touch operation and clinically tested for accuracy
- Recalls up to 99 prior readings
- World Health Organization (WHO) indicator for BP classification

Includes : Four AA batteries, DC power port and nylon storage case

Cuff Range: Large Adult: 13" - 17"*

Warranty: Two Year

**Range refers to the arm circumference that the cuff is designed to operate on*



WristMate™ Premium Digital Blood Pressure Monitor

Model : **HM-55*** (Box)

- Automatic blood pressure and pulse rate monitor that measures from the wrist
- Displays measurements in bar graph form for easy evaluation
- Features an irregular heartbeat detector
- Recalls up to 99 prior readings
- Includes World Health Organization (WHO) indicator for blood pressure classification

Includes : Two AAA batteries and hard shell storage case

Cuff Range: 5.3" - 7.5"*

Warranty: Two Year

**Range refers to the wrist circumference that the cuff is designed to operate on*



HealthMate® Premium Digital Blood Pressure Monitor

Model : **HM-35*** (Box)

- Fully automatic blood pressure and pulse rate monitor
- Displays measurements in bar graph form for easy evaluation
- Features an irregular heartbeat detector
- Recalls up to 99 prior readings
- World Health Organization (WHO) indicator for BP classification

Includes : Four AA batteries, DC power port and zippered storage case

Cuff Range: Adult: 9" - 13"*

Warranty: Two Year

**Range refers to the arm circumference that the cuff is designed to operate on*



Traditional Home Blood Pressure Set

- D-Ring in cuff enables easy self-fastening
- Black cotton cuff with owner ID and artery indicator label

Includes : Cotton cuff with attached single head stethoscope

Cuff Range: Adult : 10"-16"* , Large Adult : 13.5" - 20"*

Warranty: Lifetime Calibration

**Range refers to the arm circumference that the cuff is designed to operate on*

81 Adult Cuff (Box)

81-OB Large Adult Cuff (Box)



Temporal / Ear Thermometer

Model : **DT-28*** (Clamshell)

- Dual-Use Temporal and Ear Infrared Thermometer
- Instant one second readings
- 360° Rotating Tip for Optimal Placement
- Measures in both Fahrenheit or Celsius
- Memory of 10 Prior Readings
- Fever Alarm

Includes: CR2032 button battery
Warranty: One Year



Digital Thermometer

Model : **DT-6*** (Case)

- For use orally, rectally, or under the arm
- Measures in both Fahrenheit or Celsius
- 60 second reading
- Recalls last reading

Includes: Five probe covers, LR41 button battery
Warranty: One Year



Infrared Forehead Thermometer

Model : **DT-29*** (Box)

- Non-Contact, Fast and Accurate Readings
- Measures Body, Surface and Ambient Temperatures
- Backlit Digital LCD Display Screen
- High Temperature Alarm
- Measures in Both Fahrenheit and Celsius
- Features 32 data memory

Includes: Two AAA Batteries
Warranty: One Year



Neon Digital Thermometer

Model : **DT*** (Box)

- For use orally, rectally, or under the arm
- Measures in both Fahrenheit or Celsius
- 60 second reading
- Recalls last reading

Includes: Five probe covers, LR41 button battery
Warranty: One Year



Fingertip Pulse Oximeter

Model : **456*** (Box)

- Fingertip measurement of SpO₂% (Oxygen Saturation) and pulse rate
- Red LED display screen
- Up to 30 hours of continuous use with auto shut-off feature



BLK
Black

CGY
Cool Grey

HPK
Hot Pink

PUR
Purple

Includes: Nylon lanyard
Two AAA batteries
Soft silicone protective case
Warranty: One Year

Ice / Hot Water Bottles

- Made of natural rubber latex
- For use with cold and hot liquids
- Clear box packaging

- | | |
|-------------|-----------------------------------------|
| 4220 | 2 qt. fountain syringe kit |
| 4200 | 2 qt. Ice/Hot water bottle with stopper |
| 4210 | 2 qt. Ice/Hot stopperless water bottle |



4220



4200



4210

14" Rubber Invalid Ring

Model : **4214** (Clear Box)

- Made of natural rubber latex
- Features a metal valve
- 14" in diameter (35cm)



16" Rubber Invalid Ring

Model : **4216** (Clear Box)

- Made of natural rubber latex
- Features a metal valve
- 16" in diameter (40cm)



Standard Disposable Penlight

- Clip activated with standard illumination
- Plastic construction

200	Poly Bag
S200	Slide Pack
200-AST	Cylinder (22)

				
BLK Black	BUR Burgundy	CBL Ciel	HPK Hot Pink	HUN Hunter
				
NAV Navy	N-BLU Neon Blue	N-GRN Neon Green	N-PNK Neon Pink	PUR Purple
				
ROY Royal	TEA Teal	WHT White		

ROY
Royal



Pupil Gauge Disposable Penlight

- Pupil gauge imprint
- Clip activated
- Standard illumination

210	Poly Bag
S210	Slide Pack
210-CYL	Cylinder (22) White only
210-AST	Cylinder (22) Asst. Colors



210
White
with
Royal Print



210-HPK
Hot Pink
with
Purple Print






210-NAV
Navy
with
White Print



Pearlescent Disposable Penlight

- Clip activated
- Plastic construction with a "Pearlescent" finish
- Standard illumination

203	Poly Bag
S203	Slide Pack
203-AST	Cylinder (22)

	
AME Amethyst	RUB Ruby
	
SAP Sapphire	

RUB
Ruby



QuickLite™ Pupil Gauge Penlight

- Pupil gauge imprint
- Push button activated
- Metal construction
- Standard illumination
- Includes two AAA batteries





229	Poly Bag
S229	Slide Pack
229-CYL	Cylinder (20)



QuickLite™ Penlight

- Button activated
- Metal construction
- Standard illumination
- Includes two AAA batteries

222	Poly Bag
S222	Slide Pack
222-AST	Cylinder (20)

	
BLK Black	HPK Hot Pink
	
PUR Purple	SIL Silver

WHT
White



German Diagnostic Penlight

- Made in Germany
- Push button activated
- All-metal construction
- Standard illumination
- Includes two AAA batteries

280 Poly Bag
S280 Slide Pack



PocketLite™ Penlight

- Clip activated
- Plastic construction
- Standard illumination
- Includes two AAA batteries

260 Poly Bag
S260 Slide Pack
260-AST Cylinder (20)



Bright LED Penlight

- Button activated
- Metal construction
- Bright LED illumination*
- Includes one AAA battery


204 Poly Bag
S204 Slide Pack
204-AST Cylinder (20)



Soft LED Pupil Gauge Penlight

- Pupil gauge imprint
- Button activated with metal construction
- Soft LED illumination
- Includes one AAA battery


214 Poly Bag
S214 Slide Pack
214-AST Cylinder (20)



3-in-1 Utility Pen

- Button activated
- Bright LED Illumination*
- Easily converts from penlight to writing pen
- Features a handy stylus for use on touchscreen phones and tablets
- Includes battery

233 Poly Bag
S233 Slide Pack
233-AST Cylinder (20)



Easily flips from pen to penlight

Standard Pocket Otoscope

Model : **3420** (Slide Pack)

- Twist activated
- Dual-focused standard illumination with 2.5X magnification
- Includes five sizes of reusable specula (cold autoclavable)
- Otoscope head can be removed to convert into penlight
- Includes two AAA batteries




ROY
Royal



Mini Otolite™ Otoscope

Model : **3430** (Slide Pack)

- Clip activated
- Dual-focused standard illumination with 2X magnification
- Includes one reusable specula (cold autoclavable)
- Includes two sizes of disposable specula
- Otoscope head can be removed to convert into penlight
- Includes two AAA batteries




WHT
White



Xenon Pocket Otoscope

Model : **3470** (Slide Pack)

- Twist activated
- 6000 candle power xenon lamp with 3X magnification
- Comes with one reusable specula (cold autoclavable)
- Includes two sizes of disposable specula and coverter by Welch Allyn®
- Otoscope head can be removed to convert into penlight/flashlight
- Includes two AAA batteries



BLK
Black



Taylor Percussion Hammer

Model : **25** (Slide Pack)

- Intended to elicit patellar reflexes and myotatic reflexes
- Tapered handle tip is designed to elicit plantar reflexes
- 7½" chrome plated handle & PVC bumper



PUR
Purple



Buck Neurological Hammer

Model : **26**

- Designed to elicit myotatic reflexes and neurological responses
- Dual PVC bumpers for reflex tests
- Removable neurological pin
- Removable neurological brush
- 6¾" chrome plated handle
- Slide pack packaging



Queen Square Hammers

- Perfectly suited for eliciting a full range of reflexes
- Tapered handle base to enable plantar reflex testing
- Nylon plastic handle available in two lengths
- PVC bumper ring
- Poly bag packaging

23 14 inches
23-S 11½ inches



Babinski Telescoping Reflex Hammer

Model : **24**


- Features a 14" telescoping handle
- Enables different percussion forces
- PVC bumper ring
- Convenient to carry in pocket
- Tapered handle base to enable plantar reflex testing
- Slide pack packaging



Economy Taylor Percussion Hammer

Model : **28**

- Same features as model 25, but available in burnt orange only
- Poly bag packaging



Ergonomic Diagnostic Hammer

Model : **22**

- Combines features of the Taylor Percussion Hammer with a neurological hammer
- Retracting filament mounted at the base of the handle to elicit nervous system responses
- 6" anodized handle
- Slide pack packaging



Trömner Reflex Hammer


Model : **27-BLK**

- 3-in-1 dual mallet reflex hammer
- Larger mallet elicits tendon stretch reflex, smaller mallet elicits myotatic reflexes
- HDP plastic handle designed for a comfortable grip
- Tapered handle tip enables plantar reflex testing
- Length: 9", Weight: 5.5 oz.
- Poly bag packaging



Tuning Forks

- C-1024 intended for higher-frequency hearing tests
- C-512 & C-256 intended for Rinne and Weber hearing tests
- C-128 intended for vibration and bone fracture tests
- C-128 & C-256 weighted to assure correct frequency
- Made of aluminum



Economy EKG Calipers

Model : **36**

- Designed to measure electrocardiogram waveforms
- Plastic protective cover with molded ruler to aid with measurements
- Constructed of Zinc Alloy



German EKG Calipers

Model : **46**


- Designed to measure electrocardiogram waveforms
- Plastic protective cover with molded ruler to aid with measurements
- Constructed of Nickel-Plated Brass
- Made in Germany



Nurse Assist™ Clipboard

Model : **3309** (Box)


- Clock, timer, alarm and calculator (Includes L1131 battery)
- Made of heavy-duty ASB plastic



Goniometer


Model : 47

- Designed to measure the range of motion in joints
- Measures a 360° range in single degree increments
- Imprinted guide for gauging angles
- Conforms to ISOM standards
- Made from heavy-gauge clear plastic
- Measures in inches and centimeters



Protractor Goniometers


- Designed to measure the range of motion in joints
- Measures a 360° range in single degree increments
- Easy-read protractor
- Conforms to ISOM standards
- Made from heavy-gauge clear plastic
- Measures in inches and centimeters



CardioMeter™

Model : 31


- Ruler for interpreting EKG waveforms
- Made from heavy-gauge clear plastic



MediMeter™

Model : 49

- Pocket-size medical ruler in inches and centimeters
- Wound diameter / pupil dilation gauge
- Made from heavy-gauge clear plastic



Snellen Pocket Eye Chart

Model : 3909


- Measures visual acuity at a distance
- Dim : 18.5cm x 10cm



Jaeger Eye Chart

Model : 32

- Measures near vision
- Dim : 17.5cm x 12cm



Pocket Eye Chart

Model : 3908


- Measures basic visual acuity
- Dim : 16cm x 9cm



Pain Assessment Card

Model : **3910**

- Uses Wong-Baker Faces for pain rating scale
- Includes PQRST and COLDERRA pain assessment descriptions
- Dimensions: 18cm x 7cm



Disposable LifeMask® Face Shield

Model : **F5-24**

- Slim single use CPR face mask with 3M® Filtrete™ anti-viral filter
- Patented one-way valve with anti-viral filter
- Comes in a cylinder of 24 Face Shields



CPR LifeMask® in Nylon Keychain Bag

Model : **M2**


- Single use CPR mask with 3M Filtrete™ anti-viral filter
- LifeMask® features a patented one-way valve that prevents cross-contamination
- Comes in a nylon zippered keychain bag (three colors to choose from)



LifeMask® CPR Resuscitator

Model : **M20**

- Features a latex free air cushion that conforms to facial contours
- One-way valve with anti-viral filter to prevent cross-contamination
- Comes in hard shell storage case



CPR Resuscitator

Model : **M10**

- Features a latex free air cushion that conforms to facial contours
- One-way valve with anti-viral filter to prevent cross-contamination
- Head strap
- Comes in hard shell storage case with lanyard




Cotton Gait Belt with Metal Buckle

Model : **621**

- Used for mobilizing and transferring patients
- Made of 100% cotton with a metal buckle (58" length)



SPA
Stripes - Hot Pink



Long Cotton Gait Belt

Model : **621L-WHT**

- 100% cotton belt with a metal buckle (72" length)



Cotton Gait Belt with Plastic Buckle

Model : **620-WHT**


- Same features as model 621, but with a quick release plastic buckle



Nylon Gait Belt

- Used for mobilizing and transferring patients
- Made of heavy duty nylon (58" length)

622 Plastic Buckle
623 Metal Buckle



622-HPK



623-RAI

Medical Tourniquet

Model : **4350**

- One-hand operation with a quick release buckle
- Latex free and autoclavable
- Made in Germany



ROY -
Royal

TRO -
Tropical

RAI -
Rainbow

Lister Bandage Scissors

- Made of stainless steel
- Slide pack packaging



33
3½" Bandage



43
4½" Bandage



42
4½" Gold Plated Bandage



41
4½" with Tensionrite™ Clip



44
4½" with Large Ring



60
5½" Left-Handed Bandage



73
7¼" Bandage with Large Ring



73SR
7¼" Bandage - Serrated



53
5½" Bandage



52
5½" Gold Plated Bandage



53SR
5½" Bandage - Serrated



51
5½" with Tensionrite™ Clip



54
5½" with Large Ring

5½" Stealth Finish Lister Bandage Scissor

Model : **53-STE**

- A special stealth version of the traditional Lister bandage scissor
- Made of stainless steel
- 5½" in length
- Slide pack packaging



5½" Rainbow Finish Lister Bandage Scissor

Model : **53-RAI**

- A special rainbow version of the traditional Lister bandage scissor
- Made of stainless steel
- 5½" in length
- Slide pack packaging



Premium 5½" Lister Bandage Scissor

Model : **505**

- Made from surgical-grade 420 stainless steel
- Autoclavable (143°C / 290°F)
- 5½" in length
- Five Year Warranty
- Box packaging



Specialty / Dressing Scissors

- Made of stainless steel
- Slide pack packaging



GripMate™ 7½" Lister Bandage Scissors

Model : **94**

- Serrations in the handle are perfect for gripping or twisting
- Hook for opening injection (lidocaine) vials
- Made of stainless steel
- Slide pack packaging
- Autoclavable (143°C / 290°F)



GripMate™ 5½" Lister Bandage Scissor

Model : **84**

- Serrations in the handle are perfect for gripping or twisting
- Hook for opening injection (lidocaine) vials
- Made of stainless steel
- Slide pack packaging
- Autoclavable (143°C / 290°F)



5½" Olsen-Hegar Needle Holder / Scissor

Model : **503**

- A combination instrument that features both a needle holder and a built-in scissor
- Made of stainless steel
- 5½" in length



7½" Rochester Peaen Forceps

Model : **570**

- Features serrated jaws and a locking handle
- Made of stainless steel
- 7½" in length
- Slide pack packaging



Kelly Forceps

- Made of stainless steel
- Slide pack packaging



5½" Stealth Finish Kelly Forceps

Model : **500-STE**

- A special stealth version of the traditional Kelly forcep
- Made of stainless steel
- 5½" in length
- Slide pack packaging



5½" Rainbow Finish Kelly Forceps


Model : **500-RAI**

- A special rainbow version of the traditional Kelly forcep
- Made of stainless steel
- 5½" in length
- Slide pack packaging



Specialty Forceps & Probe

- Made of stainless steel
- Slide pack packaging



Anatomy Dissecting Kit

Model : **AK-1**

The AK-1 is our deluxe style dissection kit. It comes in a velour lined zippered carrying case (**AK-CASE**) and includes:

- 265** Two Straight Teasing Needles with Handles
- 1560** 5½" Thumb Dressing Forceps
- 1480** 4½" Sharp Splinter Forceps
- 148** 4½" Iris Scissor
- 158** 5½" Dressing Scissor (sh/sh)
- 292** Magnifying Glass
- 1495** Huber Mall Probe
- 264** Fish Hooks (3 Pack)
- 274** S/S Scalpel Handle (No. 4)
- 273** S/S Scalpel Handle (No. 3)
- 278-21** Two Scalpel Blades (No. 21)
- 278-10** Scalpel Blade (No. 10)



Standard Dissecting Kit

Model : **DK-1**

An upgraded version of the VK-1. The DK-1 features a hard plastic organizer case (**DK-CASE**) with fitted slots to protect the instruments. Kit includes:

- 265** Two Straight Teasing Needles with Handles
- 147** Student Dissecting Scissor
- 294** Glass Medicine Dropper
- 1490** 4½" Thumb Dressing Forceps
- 273** S/S Scalpel Handle (No. 3)
- 278-10** Scalpel Blade (No. 10)
- 299** Set of Six Dissecting Pins
- 295** Clear Plastic Ruler



Student Dissecting Kit

Model : **VK-1**

Ideal for dissection classes, all of the VK-1 components come in a tri-fold eyeglass-style case (**VK-CASE**). Kit includes:

- 265** Two Straight Teasing Needles with Handles
- 294** Glass Medicine Dropper
- 1560** S/S 5½" Thumb Dressing Forceps
- 147** Student Dissecting Scissor
- 273** S/S Scalpel Handle (No. 3)
- 278-10** Scalpel Blade (No. 10)
- 299** Set of Six Dissecting Pins
- 295** Clear Plastic Ruler



5½" ListerMate® Bandage Scissor

Model : **853** (Slide Pack)

- Traditional Lister design
- Exclusive molded plastic handles in an assortment of colors
- Stainless steel blades
- Autoclavable (143°C / 290°F)

Solids:




5½" ColorMate™ Kelly Forceps

Model : **504** (Slide Pack)

- Traditional Kelly forcep design
- Color-dipped handles
- Stainless steel
- Autoclavable (143°C / 290°F)

Solids:



5½" StyleMate™ Utility Scissor

Model : **871** (Slide Pack)

- Features our unique cuff prints coated on the blades
- Surgical-grade 420 stainless steel with milled shear serrations
- Autoclavable (143°C / 290°F)



CMB Blue Camouflage



CMG Green Camouflage



CMP Pink Camouflage



DAS Simple Daisies



FSH Four Square Hearts



LAB Love and Believe*



LPC Leopard Print Cream




LPG Leopard Print Grey



LVC Leaves Cream




LVG Leaves Grey



NGD Night Garden



OBK Owls Black



OCM Owls Cream



PHB Pink Hearts Black



PRB Pink Ribbons Black*




RHB Ribbons & Hearts Black*



RPK Ribbons & Hearts Pink*



RPU Ribbons & Hearts Purple*



WAB Woody Animals Black




WAC Woody Animals Cream

7½" EMT Utility Scissor

Model : **87** (Slide Pack)

- Designed to cut through clothing or bandages
- Made from surgical-grade 420 stainless steel
- Blades feature milled shear serrations
- Autoclavable (143°C / 290°F)

Solids:



7½" Premium Fluoride Scissor

Model : **607** (Box)

- Fluoride coating provides non-stick surface for cutting tape or bandages
- Surgical-grade 420 stainless steel
- Milled shear serrations
- Autoclavable (143°C / 290°F)
- Five Year Warranty

Solids:



5½" Nurse Utility Scissor

Model : **870** (Slide Pack)

- Same features as the 87 scissor but 5½" in length

Solids:



5½" Premium Fluoride Scissor

Model : **605** (Box)

- Same features as the 607 scissor but 5½" in length

Solids:



Apparel & Accessories

- 61 Scrubs
- 62 Lab Coats & Jackets
- 65 Med Sleeves™
- 66 Socks
- 68 Eyewear
- 70 Organizers
- 72 Bags & Totes
- 74 Pens
- 74 Lanyards & ID Holders
- 78 Gifts
- 79 Watches
- 83 Jewelry



Premium Five-Pocket Unisex Scrubs

Sizes : **XXS - 5X***

- Soft "Sueded" 65% Polyester, 35% Cotton Poplin blend
- High quality 5 oz. weight fabric
- V-Neck top with a tagless label
- Three pocket top with inner cell phone pocket and two pocket pant
- Bar Tack stitching for added strength
- Adjustable twill tape drawstring waist

302 Tops
401 Pants



BLK Black
CBL Ciel



GRP Grape
HUN Hunter



NAV Navy
PEW Pewter



ROY Royal
WHT White



WIN Wine



NAV
Navy



WIN
Wine

Lab Coats

Sizes : **XS - 5X***

- 65% Polyester, 35% Cotton Twill blend
- One chest pocket & two patch pockets
- 5 button front with a notched collar
- Side openings for easy access
- Two back darts & set-in back belt

5720 Womens

5710 Mens
XS not available



Women's Fashion Lab Coats

Sizes : **XXS - 3X***

- Soft "Sueded" 65% Polyester, 35% Cotton Poplin blend
- One chest pocket, two large front patch pockets and cell phone pocket
- Four-button front
- Knit cuffs & set in sleeves
- Two back darts for shaping

5821	Black
5820	White



Consultation Jackets

Sizes : **XS - 5X***

- 65% Polyester, 35% Cotton Twill blend
- One chest pocket & two patch pockets
- Three-button front with a notched collar
- Roomy interior pockets
- Set-in back belt detail

5740 Womens

5730 Mens
XS not available



Med Sleeves™ - Knitted

Model : **MS**

- Get the two shirt look without the two-shirt restrictions
- Simply slide Med Sleeves™ over your arms to provide extra warmth when you need it
- 68% cotton, 12% nylon, 17% poly and 3% spandex blend
- Over 90 styles to choose from!



HPK
Hot
Pink

SWS
Sky &
White Stripes

CYR
Cayenne
with Ruffle



WPD
White with
Pink Paws

GPR
Grey with
Pink Ruffle

SGI
Soft Grey
Fair Isle



**More Knitted styles to choose
from at www.prestigemedical.com**

Easy Fit Med Sleeves™ - Spandex

- A special Spandex version of our Med Sleeves™
- Made to fit comfortably at the bicep
- Fast drying and rated to SPF 50 for sun protection
- Stretchable Spandex / Nylon blend
- Available in two sizes

MSX Regular
MSX-S Small

More Spandex styles to choose from at www.prestigemedical.com



TQP
Trauma Queen™
on Pink

WGR
Woodland
Camouflage
Green

DPK
Digital
Camouflage
Purple

WPK
Woodland
Camouflage
Purple

Cap Case

Model : **15**

A specially designed storage case for nurse caps. Keeps caps professionally white and clean.

Cap not included.



Heartbeat EKG Nurse Socks

Model : **377-HRB**

Ankle socks made of 80% cotton, 15% nylon, and 5% spandex.

Single pairs only.



9" Compression Socks

Model : **396-WHT**

Made of 65% nylon and 35% Lycra spandex. 10 mmHg. 9" in length.

Single pairs only.
White only.



12" Compression Socks

Made of 65% nylon and 35% Lycra spandex. 10 mmHg. 12" in length. Three colors to choose from.

380 Three Pack
397 Single Pair



BLK
Black

NAV
Navy

WHT
White

12" Printed Compression Socks

Model : **386**

Made of 65% nylon and 35% Lycra spandex. Features knitted prints. 15-18 mmHg. 12" in length. Single pairs.



CMB
Blue
Camouflage

CMG
Green
Camouflage

CMP
Pink
Camouflage

NPR
Neon
Pink Ribbon*

PRP
Pink Ribbons
on Pink*

PRW
Pink Ribbons
on White*



OBK
Owls
Black

WOS
Wild Orchid
Stripes

AWS
Aqua Sea
Stripes

ICS
Island Coral
Stripes

GYS
Grey
Stripes

BKS
Black
Stripes

12" Soft Comfort Compression Socks

Model : **387**

Made of 77% polyester, 15% nylon, and 8% Lycra spandex. Features a soft finish for extra comfort. 15-20 mmHg. 12" in length. Single pairs.



WAW
Woody
Animals White

WAB
Woody
Animals Black

RHB
Ribbons &
Hearts Black*

LPC
Leopard
Print Cream

LPG
Leopard
Print Grey

Colored Full Frame Adjustable Eyewear

Model : **5400**

Anti-scratch, anti-fog, anti-static Polycarbonate lenses that meet ANSI Z87.1 standards and protect against UV rays.



Printed Full Frame Adjustable Eyewear

Model : **5420**


Anti-scratch, anti-fog, anti-static Polycarbonate lenses that meet ANSI Z87.1 standards and protect against UV rays. Adjustable temples.



Small Frame Designer Eyewear

Model : **5440**

Anti-fog, anti-static and anti-scratch Polycarbonate lenses. Meets ANSI Z87.1 standards. Features adjustable temples.



Colored Temple Eyewear

Model : **5300**

Anti-scratch, anti-fog, anti-static Polycarbonate lenses that meet ANSI Z87.1 standards and protect against UV rays.



Visitor / Student Glasses

Model : **5900**


Made in USA by Uvex and designed for general use. Features fixed side shields for improved protection.



Splash Goggles

Model : **5600**

Fog resistant, meets ANSI Z87.1 standards. Designed to protect the wearer from liquid splash by using covered vents.



Visit our website to stay updated on New Products that we introduce throughout the year.
www.prestigemedical.com

Belt Loop Organizer Kit

Features four front pockets, including a “hook-and-loop” closure coin pocket, key chain clip and matching color tabs to secure instruments. Loop on the back fits most belts. Comes with a 5½” Lister bandage scissor, disposable pupil gauge penlight and 4-color chart pen.

- 731** Full Kit
- 730** Organizer Only (Empty)

Organizer Colors



Belt Loop Organizer DX™ Kit

Model : **736**

Features four front pockets, including a “hook-and-loop” closure coin pocket, key chain clip and matching color tabs to secure instruments. Loop on the back fits most belts. Comes with a 5½” Lister bandage scissor, 5½” Kelly forceps, disposable pupil gauge penlight and 4-color chart pen.

Organizer Colors



Belt Clip Organizer Kit

Features four front pockets, a “hook-and-loop” closure coin pocket, key chain clip and matching color tabs to secure instruments. Plastic clip on the back fits any belt or waistband. Comes with a 5½” Lister bandage scissor, disposable pupil gauge penlight and 4-color chart pen.

- 733** Full Kit
- 732** Organizer Only (Empty)

Organizer Colors



5-Pocket Organizer Kit

Features four front pockets, including a “hook-and-loop” closure coin pocket, key chain clip and matching color tabs to secure instruments. Also contains a pocket on the back. Comes with a 5½” Lister bandage scissor, disposable pupil gauge penlight and 4-color chart pen.

741-WHT	Full Kit
740-WHT	Organizer Only (Empty)




5-Pocket Utility Organizer Kit

Model : **742**

A white (blank) 5-Pocket Organizer with key chain clip and “hook-and-loop” tabs to secure instruments. Comes with a 5½” Utility scissor, a color-matched disposable penlight and 4-color chart pen.

Scissor / Penlight Colors




Organizer Belt

Model : **652**

Features one tape loop, seven pockets, two “hook-and-loop” closure tabs, quick release buckle, and 50” adjustable belt. Dimensions: 8” x 8¼”.

Organizer Colors



Padded Medical Bag

Model : **753**

Specifically designed to hold medical instruments and fully lined with water-resistant nylon and protective padding. Features multiple inner pockets and zippered compartments to enable easy storing and securing of instruments. Other features include two exterior compartments, an adjustable shoulder strap, and a heavy duty handle.

Dimensions:
11" x 9" x 5".



Nurse Car-Go™ Bag

Model : **771**

The Car-Go™ utility bag is designed for the traveling or homecare nurse. Features include a travel mug holder, front mesh pouch, cell phone pouch, removable aneroid sphygmomanometer compartment and a bottom "hook-and-loop" closure strip holds the bag in place on cloth seats or floor. Multiple inner compartments hold additional equipment and personal effects.

Dimensions: 13" x 6" x 4"



Padded Nylon Holster Set

Made of heavy duty padded nylon. Five pocket design. Features a "hook-and-loop" closure front flap to secure instruments. Affixes to a belt vertically.

Set includes

- 7½" EMT utility scissor
- 5½" Kelly forceps
- 5½" thumb dressing forceps
- QuickLite™ penlight
- Buck Knife



K765 Full Set
765 Holster Only

Large Tote Bag

Model : **705**

Water-resistant lining. "Hook-and-loop" closure top, key fob and inside pocket. Dimensions: 12" x 14" x 2".



BHH
Betty
Boop®
Hearts




STG
Stick
Nurse



Fashion Tote Bag

Model : 700

Water-resistant lining. "Hook-and-loop" closure top, key fob and inside pocket. Features our fashion prints. Dim: 12" x 14" x 2".



LPC
Leopard Print
Cream



LPG
Leopard Print
Grey



LVC
Leaves
Cream




LVG
Leaves
Grey



RHB
Ribbons and
Hearts Black*



WAB
Woody
Animals Black



WAC
Woody
Animals Cream

LPG
Leopard Print
Grey



RHB
Ribbons and
Hearts Black*

Compact Carry Case

Model : 745

6" x 9" nylon zippered carry case. "Hook-and-loop" closure side pocket. Ideal for storing instruments and accessories.

WAC
Woody
Animals Cream



WOR
Wild
Orchid

Case Colors



BLK
Black



AQS
Aqua
Sea



BUR
Burgundy



CBL
Ciel



HUN
Hunter



ICO
Island
Coral



NAV
Navy



PAS
Passion



PUR
Purple



ROY
Royal




TEA
Teal




WOR
Wild
Orchid

Case Prints



LPC
Leopard Print
Cream



LPG
Leopard Print
Grey



LVC
Leaves
Cream



LVG
Leaves
Grey



PBB
Pink Ribbons
& Butterflies*




RHB
Ribbons and
Hearts Black*




RPK
Ribbons and
Hearts Pink*



RPU
Ribbons and
Hearts Purple*



WAB
Woody
Animals Black



WAC
Woody
Animals Cream

Liquid Syringe Pen

Model : **348**

Red liquid with black ink in a fun syringe pen.


Counter Display of 60 Pens



3-in-1 Utility Pen

- Writing pen that easily converts to a penlight
- Penlight is button activated with bright LED illumination*
- Stylis for use on touch screens
- Includes battery

233 Poly Bag
S233 Slide Pack
233-AST Cylinder (20)



Flips from pen to penlight



Four-Color Chart Pen

Plastic pen with black, blue, green, and red ink. Features a "no-slip" grip and pocket clip.

Plastic button to change color.

444 Poly Bag
S444 Slide Pack
444-CYL Cylinder (30)



Standard ID Holder

Model : **16**

Standard ID holder with a back bulldog clip for easy fastening to apparel or layards.



Two Way ID Holder

Model : **12**

Displays an ID either vertically or horizontally using a detachable button snap.



Magnetic ID Holder

Model : **18**

Standard ID holder with a two-piece magnetic back for easy fastening.



Basic Lanyard

Model : **6045**

Standard metal badge holder with a durable bulldog clip and a break-away clasp.



BLK
Black



PNK
Pink



PUR
Purple



PUR
Purple


Soft Brights Lanyard

Model : **6046**

Standard metal badge holder with a durable bulldog clip and a break-away clasp.



AQS
Aqua
Sea



ICO
Island
Coral



WOR
Wild
Orchid



AQS
Aqua Sea


Retractable Lanyard

Model : **6053**


Detachable Retracteze™ on a lanyard with quick release buckle. Extends up to 23". Breakaway safety clasp.



BLK
Black
Retracteze



BBW
White
Retracteze



BRB
Ruby
Retracteze



PUR
Violet
Retracteze



ROY
Black
Retracteze



PUR
Violet
Retracteze

ID Avenue™ Ribbon Lanyard with Badge Reel

Model : **BEG-0002**

- Removable badge reel with 36" retractable nylon cord
- Soft washable polyester with a breakaway safety clasp
- Hangs 19" from the back of the neck



CMB
Blue
Camouflage



CMG
Green
Camouflage



CMP
Pink
Camouflage



MNS
Monsieur
Monkey



OZO
Ozzy
Owl



ROH
Ribbons
of Hope*



SPR
Spirit of
Pink Ribbon*



WRL
Warhol



SPR
Spirit of
Pink Ribbon

Retracteze™ ID Holder

Model : S13

Retracteze™ ID badge holders extend to 23".

A metal clip on back enables easy fastening to apparel. Button-snap ID holder can also hold keys.



BLK
Black



ROY
Royal



WHT
White



VIO
Violet



RUB
Ruby



AUT
Autism



N-BLU
Neon Blue



N-GRN
Neon Green



N-PNK
Neon Pink



SVR
Silver



PRN
Pink Ribbon RN*



CBB
Love a CNA



HBB
Hearts



BFY
Butterflies



FSH
Four Square Hearts



LAB
Love and Believe*



OBK
Owls Black



OCM
Owls Cream



RPK
Ribbons & Hearts Pink*



RPU
Ribbons & Hearts Purple*



PBB
Pink Ribbons*



NGD
Night Garden




LVC
Leaves Cream



LVG
Leaves Grey



WAB
Woosy Animals Black



WAC
Woosy Animals Cream



FLN
Frog (Love a Nurse)



EKG
EKG Heart

Retracteze™ ID Holder Countertop Display

The Retracteze ID Holder Countertop display includes 60 Retracteze ID Holders (Model : 13). Individually packaged in poly bags.

DIS-13 Assorted designs (Neon excluded)

DIS-N-13 Neon colors only (Blue, Green, Pink)



ID Avenue™ Badge Reel

Model : **BEG-0001**

A spinning alligator clip back allows ID to always hang straight. 36" retractable nylon cord. Plastic snap holds any ID securely.



Deluxe Retracteze™ ID Holder

Model : **S14**

Deluxe 3-D ID badge holders extend to 23". Bulldog-style clip on back allows easy fastening to apparel. The button-snap ID holder can also hold keys.



Ceramic Graduation Lamps

Ceramic Florence Nightingale lamp with candle. Available in two sizes.

Large Dimensions : 3¼" x 6¼".

Small Dimensions : 2½" x 2½".

1890 Large

1891 Small



Battery Graduation Lamp

Model : **1897**

Large ceramic Florence Nightingale lamp with a battery operated candle.

AAA battery included. Dim: 3¾" x 6¾".



Retractable Tape Measure

Model : **45-WHT**

Measures in both inches (60") and centimeters (150cm). Easily rolls out and retracts at the push of a button.



Dr. Scrubz Bear

Model : **1841-TEA**

Measures 10" tall. Comes wearing teal scrubs with a matching scrub mask and cap.



Custom Patches

Prestige Medical specializes in custom medical patches in a wide range of styles.

Please contact customer service for details.



PROFESSIONAL WATCHES

Specially crafted for medical professionals, Prestige Medical time pieces are designed for a variety of workplace applications. All analog watches feature quartz movement with a battery included (SR626SW).

Cyber Gel™ Scrub Watch

Model : 1777

Silicone band.
Plastic case.
Plastic lens.



BLK
Black

PNK
Pink

PUR
Purple

Braided Band Scrub Watch

Model : 1889-HPK

Silicone band.
Plastic case.
Plastic lens.



SportMate™ Scrub Watch

Model : 1888

Silicone band.
Plastic case.
Plastic lens.



BLK
Black

WHT
White

PUR
Purple

N-PNK
Neon
Pink

N-BLU
Neon
Blue

Basic Scrub Watch

Model : 1770-WHT

Silicone band.
Plastic case.
Plastic lens.



Student Scrub Watch

Model : 1769

Plastic band.
Plastic case.
Plastic lens.



Neo Retro Scrub Watch

Model : 1990

Two-color plastic band.
Chrome and plastic case.
Plastic lens.



Classic Watch

Model : 1700

Leather band.
Chrome case.
Glass lens.



Easy Reader Watch

Model : 1736

Leather band.
Chrome case.
Glass lens.



Medical Symbols Watch

Model : 1799*

Leather band.
Chrome case.
Glass lens.
Available in either black or white.



Braided Band Fashion Watch

Model : 1778

Silicone band.
Chrome case.
Glass lens.
6 colors to choose from.



BLK
Black

HPK
Hot Pink

PAC
Pacific

PUR
Purple

WHT
White

Deluxe Classic Watch

Model : 1756

Leather band.
Chrome case.
Glass lens.



Two-Tone Classic Watch

Model : 1762

Leather band.
Chrome case.
Glass lens.
Available in a large or small size.



RDL
Large

RDS
Small

Nylon-Band Fashion Watch

Model : 1781

Nylon band.
Chrome case.
Glass lens.
Features medical symbols on the face.



PKS
Pink

PPW
Purple

Glow-Face Watch

Model : 1711

Leather band.
Chrome case.
Glass lens.
Features a glow-in-dark face.



High Fashion Leather Watch

Model : **1623**

Leather band.
Chrome case.
Glass lens.
Three colors
to choose
from.



Expansion Band Watch

Model : **1631**

Chrome band.
Chrome case.
Glass lens.
Features
an elastic
expansion
band to fit
any wrist.



Lanyard Watch

Model : **1789-BLK**

Black nylon lanyard with a break-
away safety clasp. Chrome case.
Glass lens.



Glow-Face Lapel Watch

Model : **1741**


Chain band. Chrome case. Glass
lens. Glow-in-the-dark face.
Features a locking pin back.



Link Lapel Watch

Model : **1740**

Linked bracelet band.
Chrome case. Glass lens.
Features a locking pin back.



Visit our website to stay updated
on New Products that we introduce
throughout the year.

www.prestigemedical.com

Insignia Tacs

Plastic clutch back.
Two per package.



90
Caduceus



91
Registered Nurse
Caduceus



92
Licensed Vocational
Nurse Caduceus



93
Licensed Practical
Nurse Caduceus



97
Certified Nursing
Assistant Caduceus



906
Medical
Assistant



910
Registered
Nurse



911
Bachelor of
Science Nursing



920
Licensed
Vocational Nurse



930
Licensed
Practical Nurse



960
Emergency
Medical Technician



975
Certified
Nursing Assistant

Professional Tacs

Plastic clutch back.
One per package.



992
Pink Ribbon*



904
Praying Hands



905
Autism



907
Eye Chart



900
Believe in
Miracles



938
Guardian
Angel



996
Tooth
Character



980
Cardio Pulmonary
Resuscitation



984
Gold Nightingale
Lamp



990
Large Gold
Caduceus

Badge Tacs

Plastic clutch back.
One per package.



9400
Pink Ribbon - Silver*



9392
Pink Ribbon*



9342
Love A Nurse



9346
We Care



9310
Registered Nurse



9326
Certified
Nursing Assistant



9330
Licensed
Practical Nurse



9314
Guardian Angel



9381
EKG Heart



9382
Registered Nurse



9390
Caduceus Blue

Magnetic Pink Ribbon Badge Tac

Model : 9412*

Two-piece magnetic
back eliminates
piercing of garment.
Crystal studded face.
One per package.



Professional Cloisonné Emblem Pins



Individually handcrafted professional cloisonné emblem pins exhibit an artisan's attention to detail. Cloisonné is an old-world metal working process that utilizes gold or silver wire to create shapes which are then filled with finely ground glass, crystals, or precious metals and fired in a kiln at extremely high temperatures. The result is exquisite beauty and luster which enhances the pride and appearance of any medical professional. All feature a locking pin back. One per package. Listed in alphabetical order.



2080
Advanced Cardiac
Life Support



2089
Anesthesia
Assistant



2011
Bachelor of
Science Nursing



1020
Caduceus



1010
Caduceus
(Filigreed Heart)



1080
Cardio Pulmonary
Resuscitation



1014
Certified
Acupuncturist



1087
Certified Emergency
Medical Technician



1073
Certified Medical
Assistant



1074
Certified Medical
Assistant



1076
Certified
Nurses Aide



1071
Certified Nursing
Assistant



1075
Certified
Nursing Assistant



2075
Certified Nursing
Assistant, Level II



1062
Certified Occupational
Therapy Assistant



1037
Certified Pharmacy
Technician



1009
Certified Registered
Nurse Practitioner



1096
Dental
Assistant



2026
EKG Technician



1090
Emergency
Medical Technician



1086
Emergency Medical
Technician



2087
Emergency
Medical Technician



1091
First
Responder



1015
Holistic
Medicine



1072
Licensed
Nursing Assistant



1023
Licensed
Practical Nurse



1013
Licensed
Practical Nurse



1033
Licensed
Practical Nurse



1053
Licensed
Practical Nurse



1022
Licensed
Vocational Nurse



1012
Licensed
Vocational Nurse



1032
Licensed
Vocational Nurse



1042
Licensed
Vocational Nurse



1008
Massage
Therapist



2020
Master of Science
Nursing



1006
Medical
Assistant



2021
Nurse
Practitioner



1007
Nursing
Assistant



1084
Paramedic



1038
Patient Care
Assistant



2024
Patient
Service Worker



1016
Pediatric Advanced
Life Support



1035
Pharmacist Rx



1046
Pharmacist Rx



2035
Pharmacy
Technician



1058
Phlebotomy
Technician



2025
Physical Therapist
Assistant



1089
Registered Dental
Hygienist



1066
Registered
Medical Assistant



1021
Registered
Nurse



1011
Registered
Nurse



1001
Registered
Nurse



1031
Registered
Nurse



1017
Registered
Nurse Practitioner



2023
Registered
Practical Nurse



1051
Registered
Veterinary Technician



1048
Respiratory
Therapist



2012
Star of Life



1093
State Tested
Nursing Assistant



2088
Surgical
Technician



2004
Usher



1049
Veterinary
Medicine



1050
Veterinary
Technician

Custom emblem pins are available, please contact customer service for details.

The Prestige Medical Website

www.prestigemedical.com

Please visit us online today. All of our diagnostic instruments can be viewed in each of their unique color offerings by clicking on the associated color swatches. We continue to enhance and improve the online shopping experience for our dealers to make it quicker and easier than ever to place purchase orders online at any time.

Like Us on Facebook
and Follow us on Twitter.



Stethoscope Warranty



Lifetime Warranty

Prestige Medical stethoscopes are warranted against defects in manufacturing during normal use by the original purchaser for the lifetime of the stethoscope. This lifetime warranty excludes the stethoscope tubing which is warranted for a period of 5 years from the date of purchase.

Any return of a Prestige Medical stethoscope must receive prior authorization by contacting Customer Service at cs@prestigemed.com or your nearest Prestige Medical dealer.

Free Lifetime Replacement Parts

Nearly all Prestige Medical stethoscopes include Free Lifetime Replacement Parts. The parts covered with this warranty program include: Eartips (pair), diaphragm, diaphragm retaining ring, non-chill ring, and bell sleeve. Free replacement parts are limited to one of each specified part per stethoscope per year. Activating this warranty is done by registering on our website at www.prestigemedical.com.

Slat Wall Display Sign

Model : **DIS-SIGN**

Features the Prestige Medical® full-color badge logo. The "Z" hook on the back allows easy installation on any slat wall. Dimensions : 30" x 14"

Please see page 23 of the 2016 Price List for more information on this item.



Do you do Hospital Shows?



Hospital Show Program

1. Please contact our sales department for a "Hospital Show Program Application Form"
2. You must have Net 30 terms to qualify and participants must be in good credit standing at the time an order is placed.
3. We'll grant net 60 terms for each purchase order placed for Hospital Shows.
4. You can have only one Hospital Show order open at any given time.
5. Product must be returned in a resalable condition. Meaning no price stickers and packaging is not damaged in any way. If the packaging has been damaged, you will be charged a 15% restocking fee.

Group Orders?

Custom Combination Kits

If you don't see a combination kit that meets your needs, simply create your own. Prestige Medical offers hundreds of possible combinations.

Build your own kit !

Simply contact
customer service
or your sales
representative.



General Return Policies

To ensure prompt handling of your returns please follow the guidelines below:

- Returns must be sent to us freight pre-paid
- No returns without prior authorization
- The name and address of your store must be enclosed inside the box
- Written reason for the return must also be enclosed inside the box

Our returns department will make a determination of whether to issue credit, replace or repair. We also reserve the right to deduct a 15% restocking charge.

3M Littmann® Stethoscopes Return Policy

If a Littmann® stethoscope or replacement part is damaged upon receipt of your order, you may return it directly to Prestige Medical®. If you have SOLD a Littmann stethoscope or replacement part to a customer (or end user), it MUST be returned by the end user directly to Littmann®.

Littmann Repairs & Refurbishing: **(800) 292-6298** • Littmann Manufacturing Defects: **(800) 228-3957**

Healthcare Apparel Sizing Guides

Prestige Scrubs™

Measurements are actual body measurements (in inches). Garment measurements will be larger to allow for ease of movement.

	XXS	XS	S	M	L	XL	2X	3X	4X	5X
CHEST	31" – 32"	32" – 34"	35" – 37"	38" – 40"	42" – 44"	46" – 48"	50" – 52"	54" – 56"	58" – 60"	62" – 64"
WAIST	23" – 24"	24" – 26"	27" – 29"	30" – 32"	34" – 36"	38" – 40"	42" – 44"	46" – 48"	50" – 52"	54" – 56"
HIP	31" – 32"	32" – 34"	35" – 37"	38" – 40"	42" – 44"	46" – 48"	50" – 52"	54" – 56"	58" – 60"	62" – 64"
Unisex inseams are approximately 31 inches										

Men's Lab Coats & Consultation Jackets

Measurements are actual body measurements (in inches). Garment measurements will be larger to allow for ease of movement.

	S	M	L	XL	2X	3X	4X	5X
CHEST	35" – 37"	38" – 40"	42" – 44"	46" – 48"	50" – 52"	54" – 56"	58" – 60"	62" – 64"
WAIST	27" – 29"	30" – 32"	34" – 36"	38" – 40"	42" – 44"	46" – 48"	50" – 52"	54" – 56"
HIP	35" – 37"	38" – 40"	42" – 44"	46" – 48"	50" – 52"	54" – 56"	58" – 60"	62" – 64"

Women's Lab Coats & Consultation Jackets

Measurements are actual body measurements (in inches). Garment measurements will be larger to allow for ease of movement.

	XS	S	M	L	XL	2X	3X	4X	5X
BUST	33" – 34"	35" – 36"	37" – 39"	40" – 43"	44" – 47"	48" – 51"	52" – 55"	57" – 59"	61" – 62"
WAIST	24" – 25"	26" – 27"	28" – 30"	31" – 34"	35" – 38"	39" – 42"	43" – 46"	48" – 50"	52" – 54"
HIP	35" – 36"	37" – 38"	39" – 41"	42" – 45"	46" – 49"	50" – 53"	54" – 57"	59" – 61"	63" – 65"

Women's Fashion Lab Coats

Measurements are actual body measurements (in inches). Garment measurements will be larger to allow for ease of movement.

	XXS	XS	S	M	L	XL	2X	3X
BUST	31" – 32"	33" – 34"	35" – 36"	37" – 39"	40" – 43"	44" – 47"	48" – 51"	52" – 55"
WAIST	24" – 25"	24" – 25"	26" – 27"	28" – 30"	31" – 34"	35" – 38"	39" – 42"	43" – 46"
HIP	33" – 34"	35" – 36"	37" – 38"	39" – 41"	42" – 45"	46" – 49"	50" – 53"	54" – 57"

For more detailed information on our full line of Healthcare Apparel, please visit us online at :

www.prestigemedical.com


Prints for Blood Pressure Cuffs and Carry Cases




CMG
Green Camouflage




CMP
Pink Camouflage



FSO
Floral Stitch Orchid



FSB
Floral Stitch Black



LPC
Leopard Print Cream



LPG
Leopard Print Grey



LVC
Leaves Cream



LVG
Leaves Grey



PBB
Pink Ribbons & Butterflies




RHB
Ribbons & Hearts Black




RPK
Ribbons & Hearts Pink



RPU
Ribbons & Hearts Purple



WAB
Woodsy Animals Black



WAC
Woodsy Animals Cream

To view virtual
product samples in
their respective colors
and prints, please
visit us online at :

www.prestigemedical.com

INDEX

Aneroid Sphygmomanometers	28–37, 39–41	ID Avenue™ Accessories (<i>BooJee™</i>)	75, 77
Aneroid Sphygmomanometer Kits	36–37, 40–41	ID Holders	74–77
Apparel	62–69	ID Tags (Stethoscope)	22
Apparel Sizing Guides	88	Jewelry	83–85
Badge Holders	74–77	Lab Coats	62–63
Badge Tacs	83	Lanyards	75
Bags (Doctors / Nurses)	72	Littmann® Stethoscopes and Products	24–27
Basics (Stethoscopes, Aneroids, Kits)	38–41	Measuring Tools	50–51
Bears (Plush Gifts)	78	Medimeter™	51
BooJee™ Accessories (<i>ID Avenue™</i>)	75, 77	Med Sleeves™	65–66
Cap Case	66	One-Hand Aneroid Sphygmomanometers	35
Carry Cases	73	Online Ordering	86
CharMED™ Accessories	23	Organizers	70–71
Clear Sound Stethoscopes	17–18	Otosopes	48
Clinical Series Combination Kits	37	Patches (Custom)	78
Clinical Series Stethoscopes	4–11	Pens	74
Clipboard	50	Penlights	46–47
Consultation Jackets	64	Pins	84–85
CPR Masks	52	Pocket Organizers	70–71
Custom Combination Kits	87	Prestige Basics (Stethoscopes, Aneroids, Kits)	38–41
Digital Blood Pressure Monitors	42–43	Print Swatches (Cuffs and Carry Cases)	89
Digital Thermometers	44	Pulse Oximeters	45
Display Sign for Slat Wall	86	Retracteze™	76–77
Disposable Stethoscope	19	Rubber Goods	45
Dissection Kits	57	Scissors and Instruments	54–59
Dual Head Stethoscopes	5–16, 42–43	Scrubs	61
Eartips	20	Single Head Stethoscopes	10, 17–19, 42–43
EKG Calipers	50	Single-Hand Aneroid Sphygmomanometers	35
EMS (EMT) Products	52–53	Socks	66–67
Electronic Stethoscopes	24	SpragueLite® Stethoscopes	14, 42
Emblem Pins	84–85	Sprague-Type Stethoscopes	10, 12–14, 42
Eye Charts	51	Stethoscopes	2–19, 24–27, 42–43
Eyewear	68–69	Stethoscope Accessories	20–23
Forceps / Hemostats	55–56	Stethoscope Parts	20
Gait Belts	53	Tacs	83
General Information and Policies	86–87	Teaching Stethoscopes	9, 13, 25, 27
Gifts	78	Thermometers (Digital)	44
Goggles / Eyewear	68–69	Tote Bags	72–73
Goniometers	51	Tourniquets	53
Graduation Lamps	78	Transfer Belts	53
Hammers (Neurological)	48–49	Tuning Forks	50
Hip Clip™ (Stethoscope Holder)	22	Veterinary Products	11
Holsters	72	Watches	79–82
Home Health Products	42–45	Warranty Information	86
Hospital Show Program	87	Web Ordering	86