





P.O. BOX 800192, Santa Clarita, CA 91380-0192
 (800) 971-9066 • (877) 971-9066 • Fax (661) 257-0688
 sales@bhsurgical.com • www.bhsurgical.com




B&H SURGICAL
 Specializing in OB/GYN and GI Products



Product Index

DuraHolder™ IPS.....	3-4
Specialty Products..... Rubber Bands, Dressings, Tapes, etc.	5-7
Sterile Aluminum Foil Products.....	8
Sterile Stockinette/Orthopedic Products..... Bias, Tubular, Impervious	9-10
Esmark Bandage Products..... Standard, Latex-Free	11-12
Sterile Elastic Bandages..... Bioline, Tensor, Coban	13-14
Sterile Cast Padding Products..... Webri™, Sof-Rol™, Specialist™, Synthetic	15-16
Cotton Products..... Applicators, Swabs, Rolls, etc.	17-18
Sterile Gauze Products..... Sponges, X-ray Detectable, Packing Gauze	19-20
Montgomery Strap Products.....	21-22
Basins, Bowls, & Labware..... Petri Dishes, Containers, Slides, etc.	23-24
Surgical Instruments & Towel Clamps..... OR Scissors, Iris, Adson, Sparta, Utility, Backhaus, etc.	25
Minor Procedure Kits & Trays..... Suture Removal, Laceration, Circumcision, Debridement, Custom Kits, etc.	26-31



The Bioseal DuraHolder™ Instrument Protection System (IPS) was designed by a leading university hospital to prevent costly damage to their delicate neuro instruments. After years of use and improvement, the system is now available to you.

You will be immediately familiar with the durable Kimguard™ fabric that we use to form each DuraHolder™ IPS.

Bioseal's unique pouch fabrication capabilities allow you to choose among one of our standard designs or customize a pouch to fit your specific requirements.

DuraHolder™ IPS is not intended to provide a sterile barrier and should be used in conjunction with an outer wrap or container system.

DuraHolder™ IPS can be used in your validated PreVac Steam Sterilizer. We recommend a minimum cycle of 270° F for 4 minutes exposure and 20 minutes dry time. As with any product you use in your sterilization system, you should validate that it works in your system. Refer to ANSI/AAMI: St79 2006 for guidance.

Kimguard is a trademark of Kimberly-Clark and is used pursuant to a licensing agreement between Bioseal and Kimberly-Clark.




DuraHolder™ IPS

DuraHolder™ IPS


DURA100/100

The DuraHolder™ 100 is designed to help you protect and organize your delicate surgical instruments.




DURA200/100

The DuraHolder™ 200 adds the convenience of notched openings and tape closure strips.



DURA300/100


The DuraHolder™ 300 has a longer flap for added protection of long instruments.



DuraHolder™ IPS


DURA050/100

The DuraHolder™ Solo 50 is a simple pocket with a closure flap. It comes in handy to control big floppy cables or tubing.




DURA505/100

The DuraHolder™ Solo 505 is a small pocket with a closure flap, designed to help protect the smallest parts.




DURA512/100


The DuraHolder™ Solo 512 is a medium pocket with a closure.



DURA515/100

The DuraHolder™ Solo 515 is a long pocket with a closure.



Ever tried to get a rope manufacturer to give you a recommendation on how to sterilize their rope? Bioseal has sterile rope, spoons, single red tops and just about any disposable class 1 medical device you can think of.

A .17566/25 - P lastic Spoon

B .17506/24 - M etal Spoon

C .16372/25 - Traction Rope 72"

D .17351/50 - Test Tube 50 m l w/screw cap

E .17308/50 - Test Tube 8 m l w/screw cap

F .17316/20 - Test Tube 16 m l w/screw cap

G .17361/50 - Test Tube 10 cc w/screw cap

H .17393/50 - G lass Tube Red Top 10 m l

I .17391/50 - Vacutainer Tube Red Top

Sterile - Latex-Free - 1/pk - 25/cs

Sterile - Latex-Free - 1/pk - 24/cs

Sterile - Latex-Free - 1/pk - 25/cs

Sterile - Latex-Free - 1/pk - 50/cs

Sterile - Latex-Free - 1/pk - 50/cs

Sterile - Latex-Free - 1/pk - 50/cs

Sterile - Latex-Free - 1/pk - 50/cs

Sterile - Latex-Free - 1/pk - 50/cs

Sterile - Latex-Free - 1/pk - 50/cs



Specialty Products



Bioseal specializes in hard to find products and specialty products that Central Supply no longer has the staff or the sterilizer space to produce. Below you will find a host of seldom thought of, but always needed, specialty implements.

Safety Pins

- Sterile or non-sterile
- Large variety of sizes to fit any need

Rubber Bands

- Sterile or non-sterile
- Latex-free
- Available in a variety of sizes

Towels, Drapes & Bags

- Sterile OR towels
- Sterile drapes
- Sterile poly or paper bags

Custom Dressings

- Sterile or non-sterile
- Made to your specification or choose from one of our special dressings

Pipe Cleaners & Brushes

- Choose from a number of standard cleaning brushes including pipe cleaners and trach brushes

Tracheostomy & Umbilical Tape

- Can be used as an economical, latex-free, alternative to tube tamers
- Available in 1/8", 1/4", 1/2" or 3/4"

Tourniquets

- Non-sterile, latex-free

Tubing

- Sterile or non-sterile
- Latex-free
- Available in a variety of sizes



Specialty Products

Safety Pins

Sterile

15001/100	#0 gold	1/pk 100pk/cs
15102/100	#1 Small	2/pk 100pk/cs
15202/100	#2 Medium	2/pk 100pk/cs
15302/100	#3 Large	2/pk 100pk/cs

Non-sterile

15000GS	#0 gold	1g/bx
15300SS	#1	1g/bx
15200SS	#2	1g/bx
15300SS	#3	1g/bx

Rubber Bands Latex-free

Sterile

14002/100	#8	2/pk 100pk/cs
14110/100	#10	2/pk 100pk/cs
14102/100	#18	2/pk 100pk/cs
14602/100	#32	2/pk 100pk/cs
14702/100	#64	2/pk 100pk/cs

Non-sterile

14000LF	#8	5g/cs
14005LF	#16	5g/cs
14100LF	#18	1g/cs
14600LF	#32	1g/cs
14200LF	#33	1g/cs
14700LF	#64	1g/cs

Towels, Drapes & Bags

1299/50	Drape Sleeve	1/pk 50pk/cs
B1706/50	Paper Bag #6	1/pk 50pk/cs
B2923/50	Clear Bag 29" x 23"	1/pk 50pk/cs
B0606/25	Resealable Bag 6" x 6"	1/pk 25pk/cs
B1215/100	Resealable Bag 12" x 15"	1/pk 100pk/cs
1050/50	Ziplock Bag 1 qt	1/pk 50pk/cs
1002/36	Alligator Bag 11" x 13"	2/pk 36pk/cs
CSR Wrapped		
KL002W/40	OR Towel Blue Cloth	2/pk 40pk/cs
KL004W/20	OR Towel Blue Cloth	4/pk 20pk/cs
KL006W/12	OR Towel Blue Cloth	6/pk 12pk/cs

Dressings

KA1124/100	Transpore Tape 1"x24"	1/pk 100pk/cs
KA1154/25	Transpore Tape 1"x54"	1/pk 25pk/cs
KA2580/50	Telfa Pad 8"x10"	1/pk 50pk/cs
KP1101/50	Maternity Pad 11"	1/pk 50pk/cs
7868/24	Combine Roll 8"x68"	1/pk 24pk/cs
7880/20	Abd Roll 8"x18"	1/pk 20pk/cs
8542/100	Orthopedic Felt 1/8"x 4" x 2.5"	1/pk 100pk/cs
11212/24	Bridal Veil 12"x 12"	1/pk 24pk/cs

Tourniquets

Non-sterile Latex-free Bulk (White)		
4918RG/125	1" x 18"	

Penrose Drain

4311/24	1" x 36"	1/pk 24pk/cs
4321/24	1/2" x 36"	1/pk 24pk/cs

Umbilical Tape

6491/72	1/4" x 90"	1/pk 72pk/cs
6812/100	1/8" x 12"	1/pk 100pk/cs
6824/120	1/8" x 24"	1/pk 120pk/cs

Trach Tape

6524/72	1/2" x 24"	2/pk 72pk/cs
6552/72	1/2" x 18"	2/pk 72pk/cs
6512/72	1/2" x 36"	2/pk 72pk/cs

Brushes

16059/24	Scrub Brush	1/pk 24pk/cs
17401/48	Toothbrush Adult	1/pk 48pk/cs
16005/50	Pipe Cleaners 6"	5/pk 50pk/cs
16011/24	Trach Brush 12" x 1/2"	1/pk 24pk/cs

Tubing

Latex-free		
NST030/20	Neuro Suction Tubing 30"	1/pk 20pk/cs
NST035/20	Neuro Suction Tubing 15"	1/pk 20pk/cs

Miscellaneous

KD202/50	Lube Jelly 3 g	2/pk 50pk/cs
----------	----------------	--------------

Cotton Roll with CSR Wrap

8460W/10	1/2 lb Roll	10/cs
8401W/05	1 lb Roll	5/cs




Sterile Aluminum Foil Products



Bioseal produces a full line of Sterile Aluminum Foil Sheets to assist with the allograft process. These sheets assist the orthopedic surgeon in structuring the bone material prior to implantation. All are packages in our durable peel pouch to ensure sterility.

17024/20	4" x 4"	1/pk	20pk/cs
17025/50	9" x 10 3/4"	1/pk	50pk/cs
17050/50	10" x 13"	2/pk	50pk/cs
17054/50	4" x 4"	1/pk	50pk/cs
17056/50	6" x 6"	1/pk	50pk/cs



As with all Bioseal products, these sterile orthopedic products are ready for use and will help you control costs and save your staff valuable time. The Bioseal Stockinette line is designed to provide a drape for orthopedic surgery.

Pre-Rolled Tubular Stockinette

- Available in single- or double-ply
- Your choice of flat fold or donut roll
- A wide variety of sizes ensures proper fit
- 100% cotton

Pre-Rolled Pull-Tab Stockinette

- Double-ply stockinette with the same features as our pre-rolled tubular stockinette
- Just pull the tabs and the stockinette easily rolls down the extremity

Pre-Rolled Bias Stockinette

- 100% cotton
- True 45 Bias cut, not cheaper cross cut

Impervious Stockinette

- A single ply of waterproof kraton plastic over polyester stockinette
- Latex-free cover

Shoulder Suspension Kit

- Mix and match Stockinette suspension sleeves of various types and lengths
- Choose your elastic bandage and other accessories

Bioseal Suture Retriever

- The knurled handles on our suture retriever have received rave reviews and acceptance from surgeons
- These suture retrievers are designed for single patient use

Custom Sizes

At Bioseal we can develop product lines to fit your specifications. If you need a size or style that you don't see here, give us a call.



Sterile Orthopedic Products

Pre-Rolled Tubular Stockinette

Sterile 1-ply Latex-free

1221/25	2" x 36"	25/cs
1231/25	3" x 36"	25/cs
1248/25	4" x 48"	25/cs
1268/25	6" x 48"	25/cs
1260/25	6" x 60"	25/cs
1282/25	8" x 72"	25/cs

Sterile 2-ply Latex-free

1321/25	2" x 36"	25/cs
1331/25	3" x 36"	25/cs
1348/25	4" x 48"	25/cs
1368/25	6" x 48"	25/cs
1360/25	6" x 60"	25/cs
1282/25	8" x 72"	25/cs

Pre-Rolled Pull-Tab Stockinette

Sterile Latex-free

1448/25	4" x 48"	25/cs
1460/25	6" x 60"	25/cs

Shoulder Suspension Kit 9" Impv

SSK119/06

Shoulder Suspension Sleeve	Qty/Kit	1
Bandage Cohesive **Latex-free**	1	1
Traction Rope	1	1
"S" Hook for Traction Kit	1	1
CSR Wrap 20" x 20"	1	1
Pouch Header 10" x 15"	1	1

Shoulder Suspension Kit

SSK149/06

Shoulder Suspension Sleeve, Impervious	Qty/Kit	1
Bandage Coban	1	1
7ft Traction Rope	1	1
"S" Hook for Traction Kit	1	1
CSR Wrap 20" x 20"	1	1
Pouch Header 10" x 15"	1	1

Bioseal Suture Retriever

RETRIEVER/5	5/cs
-------------	------

Pre-Rolled Bias Stockinette

Sterile Latex-free

1124/25	2" x 4 yds	25/cs
1134/25	3" x 4 yds	25/cs
1144/25	4" x 4 yds	25/cs
1164/25	6" x 4 yds	25/cs
1184/25	8" x 4 yds	25/cs

Sterile Latex-free

1125/25	2" x 5 yds	25/cs
1135/25	3" x 5 yds	25/cs
1145/25	4" x 5 yds	25/cs
1165/25	6" x 5 yds	25/cs
1185/25	8" x 5 yds	25/cs

Impervious Stockinette

Sterile Latex-free

1508/25	Small	6" x 30"	25/cs
1518/25	Medium	9" x 48"	25/cs
1528/25	Large	12" x 48"	25/cs
1538/25	X-Large	12" x 60"	25/cs

Shoulder Suspension Kit

SSK100/06

Shoulder Suspension Sleeve	Qty/Kit	1
Bandage Cohesive **Latex-free**	1	1
7ft Traction Rope	1	1
"S" Hook for Traction Kit	1	1
CSR Wrap 20" x 20"	1	1
Pouch Header 10" x 15"	1	1

Shoulder Suspension Kit

SSK110/06

Shoulder Suspension Sleeve, Small Impervious	Qty/Kit	1
Bandage Cohesive **Latex-free**	1	1
7ft Traction Rope	1	1
"S" Hook for Traction Kit	1	1
CSR Wrap 20" x 20"	1	1
Pouch Header 10" x 15"	1	1



Esmark Bandage Products



Bioseal Esmark Bandages are available in traditional latex or latex-free. These bandages provide strong even pressure for limb exsanguination. Bioseal sterile bandages are ready for use in the operating room. We use high quality peel pouches and state-of-the-art radiation sterilization for safe and economical sterilization.

STANDARD

Esmark Bandage

- High-quality, protein-free, latex
- No bothersome inner lining

LATEX-FREE (LF)

LF Esmark Bandage

- Latex-free alternative to our standard esmark
- No bothersome inner lining



Esmark Bandage Products

STANDARD

Esmark Peel Pouch

Sterile Standard Amber Latex

4333/20	3" x 3 yds
4334/20	3" x 4 yds
4343/20	4" x 3 yds
4344/20	4" x 4 yds
4363/20	6" x 3 yds
4364/20	6" x 4 yds

Esmark CSR Wrapped

Sterile Standard Amber Latex

4333W/20	3" x 3 yds
4334W/20	3" x 4 yds
4343W/20	4" x 3 yds
4344W/20	4" x 4 yds
4363W/20	6" x 3 yds
4364W/20	6" x 4 yds

Esmark

Non-sterile Standard Amber Latex

4333NS/20	3" x 3 yds
4334NS/20	3" x 4 yds
4343NS/20	4" x 3 yds
4344NS/20	4" x 4 yds
4363NS/20	6" x 3 yds
4364NS/20	6" x 4 yds

LATEX-FREE

LF Esmark Peel Pouch

Sterile Light Blue

20/cs	4933/20	3" x 3 yds	20/cs
20/cs	4934/20	3" x 4 yds	20/cs
20/cs	4943/20	4" x 3 yds	20/cs
20/cs	4944/20	4" x 4 yds	20/cs
20/cs	4963/20	6" x 3 yds	20/cs
20/cs	4964/20	6" x 4 yds	20/cs

LF Esmark CSR Wrapped

Sterile Light Blue

20/cs	4933W/20	3" x 3 yds	20/cs
20/cs	4934W/20	3" x 4 yds	20/cs
20/cs	4943W/20	4" x 3 yds	20/cs
20/cs	4944W/20	4" x 4 yds	20/cs
20/cs	4963W/20	6" x 3 yds	20/cs
20/cs	4964W/20	6" x 4 yds	20/cs

LF Esmark

Non-sterile Light Blue

20/cs	4933NS/20	3" x 3 yds	20/cs
20/cs	4934NS/20	3" x 4 yds	20/cs
20/cs	4943NS/20	4" x 3 yds	20/cs
20/cs	4944NS/20	4" x 4 yds	20/cs
20/cs	4963NS/20	6" x 3 yds	20/cs
20/cs	4964NS/20	6" x 4 yds	20/cs



Sterile Elastic Bandage Products



Bioseal offers a wide variety of elastic bandages both standard and latex-free. All are constructed of the highest quality materials, packaged in our durable peel pouch to ensure sterility.

STANDARD

Coban®

- Adheres to itself, no need for tape or clips
- 3M "Brand Name" quality

Tensor® Premium Elastic Bandage

- Premium quality cotton and high-count rubber
- Pre-attached closure clips to ease application
- Kendall "Brand Name" quality

Bioline Standard Elastic Bandage

- High quality at an economical price

Bioline Velclose Elastic Bandage

- High-quality elastic bandage with velcro closures instead of clips

LATEX-FREE (LF)

LF Cohesive Bandage

- Adheres to itself, no need for tape or clips.
- High-quality latex-free alternative to Coban®

LF Elastic Bandage

- Latex-free at an economical price
- Cotton fiber and Kraton elastic for comfort

LF Bioline Velclose Elastic Bandage

- High-quality elastic bandage with velcro closures instead of clips

Coban is a registered trademark of 3M

Tensor is a registered trademark of The Kendall Co.



Sterile Elastic Bandage Products

STANDARD

Coban®

Sterile

4015/30	1" x 5 yds	30/cs
4025/24	2" x 5 yds	24/cs
4035/24	3" x 5 yds	24/cs
4045/24	4" x 5 yds	24/cs
4065/24	6" x 5 yds	24/cs

Tensor® Premium Elastic Bandage

Sterile Premium Grade Pre-attached Clips

4125/24	2" x 5 yds	24/cs
4135/24	3" x 5 yds	24/cs
4145/24	4" x 5 yds	24/cs
4165/24	6" x 5 yds	24/cs

Bioline Standard Elastic Bandage

Sterile Standard Grade with Clips

4225/25	2" x 5 yds	24/cs
4235/25	3" x 5 yds	24/cs
4245/25	4" x 5 yds	24/cs
4265/25	6" x 5 yds	24/cs
4260/25	6"x 10 yds	24/cs

Bioline Velcrose Elastic Bandage

Sterile with Velcro Closure

4525/25	2" x 5 yds	24/cs
4535/25	3" x 5 yds	24/cs
4545/25	4" x 5 yds	24/cs
4565/25	6" x 5 yds	24/cs

LATEX-FREE

LF Cohesive Bandage

Sterile

4026/24	2" x 5 yds	24/cs
4036/24	3" x 5 yds	24/cs
4046/24	4" x 5 yds	24/cs
4066/24	6" x 5 yds	24/cs

LF Elastic Bandage

Sterile Standard Grade with Clips

4725/25	2" x 5 yds	25/cs
4735/25	3" x 5 yds	25/cs
4745/25	4" x 5 yds	25/cs
4760/25	6" x 11 yds	25/cs
4765/25	6" x 5 yds	25/cs

LF Bioline Velcrose Elastic Bandage

Sterile Latex-free with Velcro Closure

4825/25	2" x 5 yds	25/cs
4835/25	3" x 5 yds	25/cs
4845/25	4" x 5 yds	25/cs
4865/25	6" x 5 yds	25/cs



Sterile Cast Padding Products



Bioseal produces a full line of Sterile Cast Padding to fill the needs of the orthopaedic surgeon. This includes BSN's Sof-Rol™ and Specialist™, Kendall's Webril™ line, and our own Bi-Web Cotton Cast Padding and Bioline Synthetic Cast Padding. All of our Cast Paddings are radiation sterilized to insure maximum patient safety and no risk of ETO residual.

Bi-Web

- Best value
- 100% cotton fiber for comfort
- Sterile or non-sterile

Webril™

- Sized finish for added strength
- 100% cotton fiber for comfort
- Kendall "Brand Name" quality

Bioline Synthetic

- Designed for use under new hi-tech casts

Sof-Rol™

- 100% rayon fiber
- "Brand Name" quality

Specialist™

- Pleated design for added conformability
- 100% cotton fiber for comfort
- "Brand Name" quality

Webril is a trademark of The Kendall Company
Specialist and Sof-Rol are trademarks of BSN



Sterile Cast Padding Products

Bi-Web

Sterile

3024/24	2" x 4 yds
3034/24	3" x 4 yds
3044/24	4" x 4 yds
3064/24	6" x 4 yds

Webril™

Sterile

3024k/24	2" x 4 yds
3034k/24	3" x 4 yds
3044k/24	4" x 4 yds
3064k/24	6" x 4 yds

Synthetic

Sterile

3534/20	3" x 4 yds	20/cs
3544/20	4" x 4 yds	20/cs
3564/20	6" x 4 yds	20/cs

Sof-Rol™

Sterile

3124/24	2" x 4 yds	24/cs
3134/24	3" x 4 yds	24/cs
3144/24	4" x 4 yds	24/cs
3164/24	6" x 4 yds	24/cs

Specialist™

Sterile

3224/24	2" x 4 yds	24/cs
3234/24	3" x 4 yds	24/cs
3244/24	4" x 4 yds	24/cs
3264/24	6" x 4 yds	24/cs



Bioseal saves you and your staff time and money by having a line of convenient sterile cotton goods ready for use in the operating room.

Absorbent Balls

- Available in cotton or more economical rayon
- Wide variety of sizes
- Some are even X-ray detectable

Cotton Swabs

- Large tip 8" OB swabs
- Large tip procto swabs

Cotton Tip Applicators

- 3" or 6" length available
- Packaged in 5, 10, 12, 18, 20, and 50 per pk

Cotton Rolls

- Ready for use in the OR
- 1/2 or 1 pound rolls

Cotton Squares

- Pre-cut surgical cotton
- Use for padding or dressings

Orthopedic Felt

- Pre-packaged and sterile
- Available in a variety of sizes



Cotton Products

Absorbent Cotton Balls

Sterile

8102/100	Medium	2/pk	100pk/cs
8101/100	Medium	1/pk	100pk/cs
8102/50	Medium	2/pk	50pk/cs
8103/50	Medium	3/pk	50pk/cs
8105/100	Medium	5/pk	100pk/cs
8110/100	Medium	10/pk	100pk/cs
8201/100	Large	1/pk	100pk/cs
8201/50	Large	1/pk	50pk/cs
8205/100	Large	5/pk	100pk/cs
8210/100	Large	10/pk	100pk/cs
8305/50	X-Large	5/pk	50pk/cs
8310/100	X-Large	10/pk	100pk/cs

X-Ray Detectable Absorbent Cotton Balls

Sterile

8105U/50	Small w/String	5/pk	50pk/cs
8205U/50	Large w/String	5/pk	50pk/cs
8285/50	Large 1 1/4" wide	5/pk	50pk/cs

Absorbent Rayon Balls

10120/50	Medium	20/pk	50pk/cs
10105/100	Medium	5/pk	100pk/cs
10110/100	Medium	10/pk	100pk/cs
10205/100	Large	5/pk	100pk/cs

Cotton Swabs

9702/50	Cotton Swabs 8	2/pk	50pk/cs
9712/100	OB/Gyn Swabs 8	2/pk	100pk/cs
9703/50	Procto Swabs 13 1/2"	3/pk	50pk/cs
9706/100	Procto Swabs 16 1/2"	6/pk	100pk/cs
9713/50	Procto Swabs 13 1/2"	2/pk	50pk/cs
9716/50	Procto Swabs 16 1/2"	2/pk	50pk/cs

3" Cotton Tip Applicators

9305/50		5/pk	50pk/cs
9310/100		10/pk	100pk/cs
9312M/50		12/pk	50pk/cs
9318/100		18/pk	100pk/cs
9320/100		20/pk	100pk/cs
9320P/100	Plastic Shaft	20/pk	100pk/cs
9325/100		25/pk	100pk/cs
9350/100		50/pk	100pk/cs
9399/50		100/pk	50pk/cs

6" Cotton Tip Applicators

9612/100		12/pk	100pk/cs
9620/100		20/pk	100pk/cs
9602/100		2/pk	100pk/cs
9605/100		5/pk	100pk/cs
9606L/100	Large Tip	6/pk	100pk/cs
9610/100		10/pk	100pk/cs
9625/100		25/pk	100pk/cs
9650/50		50/pk	50pk/cs

Cotton Rolls

8460W/10	1/2 lb Cotton Roll	1/pk	10pk/cs
8401/05	1 lb Cotton Roll	1/pk	5pk/cs
8468/20	Dr. Sheridan's Cotton 6" x 8" 2-ply	1/pk	20pk/cs
8401WN/05	1 lb Cotton Roll w/o Inner Liner	1/pk	5pk/cs

CSR Wrapped

8465W/10	6" x 5 ft	1/pk	10pk/cs
8402W/10	12" x 2 yds	1/pk	10pk/cs
8401W/05	1 lb Cotton Roll	1/pk	5pk/cs
8430/06	Cotton Roll 30	1/pk	6pk/cs

Cotton Squares

8421/50	2" x 2"	1/pk	50pk/cs
8441/50	4" x 4"	1/pk	50pk/cs
8413/50	1" x 3"	1/pk	50pk/cs



Bioseal offers a line of specialty gauze products designed to serve the needs of operating room staff. They are made of high-quality cotton gauze and are pre-counted to help maintain accountability. We then utilize state-of-the-art radiation sterilization for safe and economical sterilization.

Sterile Packing Gauze

- Available in 1" or 2" widths
- X-ray detectable

Sterile Gauze Sponges

- Pre-counted to insure proper accountability
- Available in 5- or 10-count packs
- 8-ply, 12-ply, 16-ply, Standard or X-ray

Sterile Tonsil Sponges

- X-ray detectable for easy retrieval
- Pre-counted in 5/pk
- Cotton filled gauze for absorbency

Sterile Stick Sponges

- X-ray detectable for easy retrieval
- Pre-counted in 5/pk to insure accountability
- For use in deep stick sponging or prepping

Sterile Dissecting Sponges

- X-ray opaque for easy retrieval
- In foam holders or pinned, Peanut, Cherry and Kittner

Sterile Dental Rolls

- Available both with or without string
- Packed 2/pk, 5/pk
- Strung are X-ray detectable

Sterile Stretch Gauze

- Available in a variety of sizes
- Conforms easily to body contours
- High cotton content for excellent absorbency

Sterile Roller Gauze

- 1", 2", 3", 4", and 6" widths
- 44 x 36 or 24 x 48 mesh
- Full rolls or cut to length



Sterile Gauze Products

Sterile Packing Gauze

X-ray Detectable

2103/25	1" x 3 yds	25/cs
2104/25	1" x 4 yds	25/cs
2105/25	1" x 5 yds	25/cs
2201/25	2" x 1 yd	25/cs
2202/25	2" x 2 yds	25/cs
2203/25	2" x 3 yds	25/cs
2204/25	2" x 4 yds	25/cs
2205/25	2" x 5 yds	25/cs

Sterile Packing Gauze

X-ray Detectable with Loop

2411L/25	4" x 1 yd	25/cs
2412L/25	4" x 2 yds	25/cs
2415L/25	4" x 5 yds	25/cs

Sterile Gauze Sponges

Pre-Counted, Standard 12-ply

7125/50	2" x 2"	5/pk	50/cs
7120/50	2" x 2"	10/pk	50/cs
7135/50	3" x 3"	5/pk	50/cs
7130/50	3" x 3"	10/pk	50/cs
7145/50	4" x 4"	5/pk	50/cs
7140/50	4" x 4"	10/pk	50/cs
7185/50	8" x 4"	5/pk	50/cs
7180/50	8" x 4"	10/pk	50/cs

Pre-Counted, X-ray Detectable 12-ply

7325/50	2" x 2"	5/pk	
7320/50	2" x 2"	10/pk	
7335/50	3" x 3"	5/pk	
7330/50	3" x 3"	10/pk	
7345/50	4" x 4"	5/pk	
7340/50	4" x 4"	10/pk	
7385/50	8" x 4"	5/pk	
7380/50	8" x 4"	10/pk	

Sterile Stretch Gauze

4415/24	1" x 5 yds	24/cs
4425/24	2" x 5 yds	24/cs
4435/24	3" x 5 yds	24/cs
4445/24	4" x 5 yds	24/cs
4465/24	6" x 5 yds	24/cs

Sterile Roller Gauze

24 x 48 Mesh

7501/25	1" x 10 yds	25/cs
7502/25	2" x 10 yds	25/cs
7503/25	3" x 10 yds	25/cs
7504/25	4" x 10 yds	25/cs
7506/25	6" x 10 yds	25/cs

44 x 36 Mesh

7601/25	1" x 10 yds	25/cs
7602/25	2" x 10 yds	25/cs
7603/25	3" x 10 yds	25/cs
7604/25	4" x 10 yds	25/cs
7606/25	6" x 10 yds	25/cs

Tonsil Sponge

Sterile X-ray Detectable, 5 On a Card

13105/50	Small	50/cs
13205/50	Medium	50/cs
13305/50	Large	50/cs
13405/50	X-Large	50/cs

Stick Sponge

Sterile X-ray Detectable, 5 On a Pin

12105/50	Small	50/cs
12205/50	Medium	50/cs
12505/50	Large	50/cs
12305/50	X-Large	50/cs

Dissecting Sponges

Sterile X-ray Detectable, 5 In a Holder

11005/50	Peanut	50/cs
15605/50	Kittner	50/cs
13505/50	Cherry	50/cs

Dental Rolls

Sterile with String

9901/50	1/pk	24/cs
9942/50	2/pk	24/cs
9905/50	5/pk	24/cs

Sterile without String

9916/50	6/pk	24/cs
9918/50	18/pk	24/cs



Bioseal Montgomery Straps are now available in easy access dispenser boxes and can be cut to any desired width, or you can choose one of our sterile pre-cut kits that include straps and twill ties.

Features & Benefits

- Direct replacement for the discontinued Dermicel® Montgomery Strap from J & J 5129
- Latex-free
- Reinforced tie holes
- Quick dressing changes save time
- No need to remove and reapply tape during dressing changes which minimizes skin trauma
- Available in dispenser boxes for cut-to-fit or convenient ready-to-apply sterile kits

Dermicel is a registered trademark of Johnson & Johnson



Montgomery Straps

Sterile Montgomery Strap Kit

5001/48	1 Hole
5002/48	2 Hole
5003/24	3 Hole
5004/24	4 Hole

Non-sterile Montgomery Strap Kit

5001NS/48	1 Hole
5002NS/48	2 Hole
5003NS/24	3 Hole
5004NS/24	4 Hole

Non-sterile Montgomery Strap Dispenser

5129BX	Box	24 Sheets
5129CS	Case	6 Boxes



Bioseal's labware products are designed to serve the needs of operating room staff. These items are conveniently packaged and ready for use. As with all Bioseal products, they will help you control costs and save your staff valuable time

Sterile Bowls, Basins & Pitchers

- 32 oz solution bowl
- 700 cc emesis basin
- 16 oz solution bowl
- 7 qt bowl
- 1200 cc pitcher grad w/handle

Sterile Specimen Containers & Cups

- 1 oz medicine cup
- 2 oz medicine cup
- 5 oz specimen cup w/lid
- 8 oz specimen cup w/lid

Sterile Tubes

- 10 cc test tube w/screw cap
- Test tube w/screw cap
- 50 ml test tube w/cap

Sterile Conectors

- 5-in-1 connector
- Simms connector
- 6-in-1 Y connector
- 6-in-1 connector

Sterile Medicine Droppers

- 1 ml plastic medicine dropper
- 1 ml glass medicine dropper

Sterile Microscope Slides

- Microscope slides
- Microscope slides w/mailer

Sterile Petri Dishes

- 100 mm x 15 mm Petri dish



Basins, Bowls, & Labware

Sterile Bowls, Basins & Pitchers

KB007W/20	Emesis Basin 700 cc	1/pk	20pk/cs
KB016W/20	Solution Bowl 16 oz CSR Wrapped	1/pk	20pk/cs
KB032W/20	Solution Bowl 32 oz CSR Wrapped	1/pk	20pk/cs
KB017W/10	Bowl 7 Quart CSR Wrapped	1/pk	10pk/cs
17511/20	Pitcher Grad w/Handle 1200cc CSR Wrapped	1/pk	20pk/cs

Sterile Specimen Containers & Cups

17521/50	Medicine Cup 1 oz	1/pk	50pk/cs
17521/100	Medicine Cup 1 oz	1/pk	100pk/cs
17521T/50	Medicine Cup 1 oz w/Toothpick	1set/pk	50pk/cs
17522/50	Medicine Cup 2 oz	1/pk	50pk/cs
17522/100	Medicine Cup 2 oz	1/pk	100pk/cs
17523/100	Medicine Cup 2 oz w/Lid	1/pk	100pk/cs
17531/50	Specimen Cup 5 oz w/Lid	1/pk	50pk/cs
KB008W/50	Solution Cup 8 oz w/Lid CSR Wrapped	1/pk	50pk/cs

Sterile Petri Dishes

17651/100	Petri Dish 100 mm x 15 mm	1/pk	100pk/cs
17240/48	Petri Dish 100 mm x 15 mm	1/pk	48pk/cs

Sterile Connectors

17601/50	Connector 5-in-1	1/pk	50pk/cs
17611/50	Connector Simms	1/pk	50pk/cs
17641/50	Connector 6-in-1 Y	1/pk	50pk/cs
17602/50	Connector 6-in-1	1/pk	50pk/cs

Sterile Medicine Droppers


17331/25	Medicine Dropper Plastic 1 ml	1/pk	25pk/cs
17340/50	Medicine Dropper Glass 1 ml	1/pk	50pk/cs

Sterile Microscope Slides

17002/50	Microscope Slides 3" x 1"	2/pk	50pk/cs
17022/50	Microscope Slides w/Mailer	1set/pk	50pk/cs

Sterile Tubes

17361/50	Test Tube 10 cc w/Screw Cap	1/pk	50pk/cs
17361/100	Test Tube 10 cc w/Screw Cap	1/pk	100pk/cs
17351/50	Test Tube 50 ml w/Cap	1/pk	50pk/cs



Bioseal offers a comprehensive line of sterile single-use surgical instruments packaged in our durable peel pouch. Other sizes and patterns available upon request.

Sterile Scissors

50/case		Length
KI011/50	OR Sharp/Blunt	5 1/2"
KI031/50	Straight Iris	4 1/2"
KI032/50	Curved Iris	4 1/2"
KI034/50	Wire Cutting	4 3/4"

Sterile Hemostatic Forceps

50/case		Length
KI081/50	Curved Mosquito	5"
KI082/50	Straight Mosquito	5"
KI084/50	Curved Kelly	5 1/2"
KI083/50	Straight Kelly	5 1/2"
KI085/50	Straight Rochester Pean	6 1/4"
KI086/50	Curved Rochester Pean	7 1/4"

Sterile Thumb Dressing Forceps

50/case		Length
KI071/50	Adson Tissue w/1x2 Teeth	4 3/4"
KI072/50	Adson Dressing, serrated	4 3/4"
KI075/50	Standard Dressing, serrated	4 3/4"
KI076/50	Standard Tissue w/1x2 Teeth	4 3/4"



Bioseal offers a wide selection of towel clamps. These clamps are available individually packaged sterile or non-sterile.

Sterile Towel Clamps

1/pk 50/case	
SP55/50	Sparta Clamp
KI105/50	Plastic Towel Clamp
KI051/50	Backhaus Towel Clamp Sm
KI055/50	Backhaus Towel Clamp Lg

Non-sterile Towel Clamps

Each	
SP55	Sparta Clamp
KI105	Plastic Towel Clamp
KI051	Backhaus Towel Clamp Sm
KI055	Backhaus Towel Clamp Lg

Sterile Needle Holders

50/case		Length
KI061/50	Webster, Smooth Jaw	5"
KI061S/50	Webster, Satin, Smooth Jaw	5"
KI062/50	Baumgartner, Serrated Jaw	5"
KI065/50	Mayo Hegar	6"



Minor Procedure Kits & Trays



BREAST BIOPSY
CHEST TUBE
CIRCUMCISION
DEBRIDEMENT
GENERAL PURPOSE
HEAD DRESSING
LACERATION
NASAL
PELVIC EXAM
PREP TRAY
SUTURE REMOVAL
SHOULDER SUSPENSION
VASECTOMY



Sigmoidoscope Kits



Bioseal's disposable sigmoidoscope and anoscope kits help eliminate the risk of cross contamination. They are 100% latex free, require no filter and eliminate the expense associated with the new cleaning requirements. Choose one of our premade kits or customize your own.

	KI510/10	KI511/10	KI512/10	KI513/10	KI514/10	KI515/10	KI516/10
W-A Sigmoidoscope #53130	•	•	•	•	•	•	•
Bulb and Bladder	•	•	•	•	•	•	•
Procto Swabs		2		4	2		2
Suction Sigmoidoscope Conmed (#0033050)				•		•	
Suction Sigmoidoscope Busse (#165)			•				
W-A Anoscope #53110					•		





Minor Procedure Kits & Trays



Breast Biopsy Tray

BBT003/20

Towel Absorbent 2-ply White
Gauze Sponge 4" x 4" 8-ply
Cotton Tip Applicators 6"
Steri Strip 1/4" x 3"
Needle 18 g x 1½"
Needle 25 g x 5/8"
Syringe 10 cc Luer Lock
Syringe 5 cc Luer Lock
Syringe 20 cc Luer Lock
Scalpel Disposable #11
Forceps Sponge 7 1/4" Plastic
Instrument Pack Pouch
Tray Small 4.74" x 5.40 x 1"
Tray 9" X 5" x 2"
CSR Wrap 24" x 24"

Qty/Kit

1
10
2
1
1
1
1
1
1
1
1
1
1
1
1



Pelvic Exam Tray Medium

PLV009/20

Cotton Ball - Large
Cotton Swab 8"
Gauze Sponge 4x4 8-ply
Glove Medium **Latex-free**
Lube Jelly Packet 3 g
Pouch Header 10" x 15"
Speculum Fiber Optic
Sterilization Indicator Blank
Tab Label
Towel Absorbent 2-ply White
Towel Poly Back - 13" x 19"
Tray 9" x 5" x 2"
CSR Wrap 24" x 24"

Qty/Kit

8
3
2
2
2
1
1
1
1
1
1
1
1



Debridement Set

DBR005/50

Adson Tissue Forceps 4¾"
Gauze Sponge Std 3" x 3" 12-ply
Pouch 6" x 10"
Tenotomy Scissors Cvd 4.5"
CSR Wrap 12" x 12"

Qty/Kit

1
6
1
1
1



Shoulder Suspension Kit 9" Impv

SSK119/06

Shoulder Suspension Sleeve
Bandage Cohesive **Latex-free**
Traction Rope
"S" Hook for Traction Kit
CSR Wrap 20" x 20"
Pouch Header 10" x 15"

Qty/Kit

1
1
1
1
1
1




Minor Procedure Kits & Trays



Circumcision Tray

CIR001/20

Adson Tissue Forceps 4 $\frac{3}{4}$ "	1
Drape Polylined 18" x 26"	1
Gauze Sponge 4" x 4" 8-ply	4
Instrument Pack Pouch	2
Iris Scissors Curved 4 $\frac{1}{2}$ "	1
Mosquito Forceps Curved 5"	2
Mosquito Forceps Str. 5"	2
OR Scissors 5 1/2" Str. S/B	1
Pouch Header 10" x 15"	1
Probe w/Eye 5"	1
Safety Pin #1 Small Steel	1
Scalpel Disposable #10	1
Sterilization Indicator Blank	1
Tab Label	1
Towel Absorbent 2-ply White	2
Tray 9" x 5" x 2"	1
CSR Wrap 24" x 24"	1



Vasectomy Tray

VST405/20

Adson Tissue Forceps 4 $\frac{3}{4}$ "	1
Allis Tissue Forceps "Baby"	2
Backhaus Towel Clamps 3 $\frac{1}{2}$ "	2
Gauze Sponge 4" x 4" 8-ply	4
Gauze Sponge Std 4" x 4" 12ply	4
Instrument Pack Pouch	3
Littauer Stitch Scissors 4 $\frac{1}{2}$ "	1
Medicine Cup 2 oz	1
Mosquito Forceps Curved 5"	2
Mosquito Forceps Str 5"	1
Needle Holder Webster	1
Pouch Header 10" x 15"	1
Scalpel Disposable #15	1
Tray 7 1/2" x 5" x 1 $\frac{1}{4}$ "	1
Paper Wrap 20" X 20"	1



Head Dressing Tray

HDT200/20

Cotton Roll	1
Kerlix Roll	2
Non Woven Sponges 4" x 4" 4-ply	10
Pouch Header 10" x 15"	1
Telfa Pad 8" x 10"	1
Wrap CSR 20" x 20"	1



Suture Removal Set

SUT302/50

Adson Dressing Forceps 4 $\frac{3}{4}$ "	1
Gauze Sponge Std 2" x 2" 8-ply	2
Instrument Pack Pouch	1
Iris Scissors Straight 4 $\frac{1}{2}$ "	1
Pouch 4" x 9"	1



General Purpose Tray

GPT006/50

	Qty/Kit
Adson Dressing Forceps 4 ³ / ₄ "	1
Alcohol Prep Pad	1
Cotton Tip Applicators 3"	2
Gauze Sponge Std 4" x 4" 12-ply	2
Instrument Pack Pouch	1
Kelly Forceps Straight 5 ¹ / ₂ "	1
Label **Expiration Date**	1
OR Scissors 5 ¹ / ₂ " Str. SB	1
Pouch 8" x 12	1
PVP Iodine Prep Pad Medium	1
Safety Pin #3 Large Steel	1
Tray 7 1/2" x 5" x 1 1/4"	1
CSR Wrap 20" x 20"	1



Instrument Set

INST03/20

	Qty/Kit
Adson Dressing Forceps 4 ³ / ₄ "	1
Alcohol Prep Pad	1
Cotton Tip Applicators 3"	2
Gauze Sponge Std 3" x 3" 12-ply	2
Instrument Pack Pouch	1
Iris Scissors Straight 4 ¹ / ₂ "	1
Label **Expiration Date**	1
Mosquito Forceps Curved 5"	1
Pouch 4" x 9"	1
PVP Iodine Prep Pad Medium	1
Safety Pin #3 Large Steel	1



Nasal Set

NSL001/20

	Qty/Kit
Dental Rolls Cotton NO Strings	2
Frazier Suction Vented 12 FR	1
Gauze Sponge Std 2" x 2" 8-ply	5
Jansen Ear Forceps Serrated	1
Nasal Speculum Plastic Adult	1
Pouch 6" x 12"	1
Wrap CSR 20" X 20"	1



Prep Tray

PRP008/20

	Qty/Kit
Pouch Header 10" x 15"	1
Towel Absorbent 2-ply White	2
Towel Poly Back - 13" x 19"	2
Winged Prep Sponges	4
Glove Medium **Latex-free**	2
Glove Paper to make Walledd	1
Tray 2 Compartments	1
Wrap CSR 24" x 24"	1



	2
10	1
	1
	1
	1
	1
	1
	1
	1
	1
	1
	1
	3
	1
	1
	1

1
1
10
5
1
1
2
1
1
1
1
1
1
1
1
2
1
1

Bioseal features a wide variety of manufacturers, including PDI, Triad, Smith & Nephew, BD, Aplicare, J & J, Kendall, 3M, Smith & Nephew, and Venigard.



Distributed by



P.O. Box 800192, Santa Clarita, CA 91380-0192 • (800) 971-9066 • (877) 971-9066 • Fax (661) 257-0688
sales@bhsurgical.com • www.bhsurgical.com