

Contents	Page
Biopsy Forceps - Single Use	
Single Use Biopsy Forceps – Cold	1
Single Use Bronchoscope Biopsy Forceps	1
Single Use Biopsy Forceps – Hot	1
Biopsy Forceps - Reusable	
Reusable Biopsy Forceps - with & without Spike	2
Reusable Bronchoscope Biopsy Forceps - without Spike	2
Reusable Biopsy Forceps - Coated Sheath	2
Hot Reusable Forceps	2
Retrieval Devices - Single Use & Reusable	
Tripod Graspers	3
Foreign Body Baskets - Elliptical & Hexagonal	3
Coagulation Probes	
Bipolar Coagulation Probes	4
Polyp Catches	
Two-Chamber Polyp Catches	4
Sclerotherapy / Injection Needles	
Injection / Sclerotherapy Needles - Single Use	5
Transbronchial Aspiration Needles	
Transbronchial Aspiration Needles - Single Use	5
Polypectomy Snares	
Polypectomy Snares - Single Use	6
Polypectomy Snares - Reusable	7
Channel Cleaning Brushes	
Reusable Channel Cleaning Brushes	8
Single Use Channel Cleaning Brushes	8
Cytology Brushes	
Cytology Brushes	8
Patient Accessories	
Privacy Pants - Single Use	9
Bite Blocks - Single Use	9
Terms & Conditions	
Terms & Conditions	10

Single Use Biopsy Forceps

»» Single-Use Convenience
»» Coated Metal Spiral Sheath

»» Stainless Steel Jaws
»» Sterile Packed

Cold Biopsy Forceps

	Description	Catalog No.	Length (cm)	Outer Diameter (mm)
	w/o spike, fenestrated ellipsoid cups	XCBF 40615	160	2.4
		XCBF 40622	230	2.4
	w/spike, fenestrated ellipsoid cups	XCBF 40715	160	2.4
		XCBF 40722	230	2.4

10/box

Bronchoscope Biopsy Forceps

	Description	Catalog No.	Length (cm)	Outer Diameter (mm)
	w/o spike, fenestrated ellipsoid cups, non-coated metal spiral	XCBF 40115	160	1.7

10/box

Hot Biopsy Forceps

	Description	Catalog No.	Length (cm)	Outer Diameter (mm)
	w/o spike, fenestrated ellipsoid cups	XHBF 40615	160	2.4
		XHBF 40622	230	2.4

10/box

Reusable Biopsy Forceps

- »» Autoclavable
- »» Durable Handle
- »» Smooth, Tapered Distal Spiral
- »» Excellent Value
- »» Hardened Stainless Steel Jaws
- »» Easy to Clean & Maintain

Cold Biopsy Forceps

	Description	Catalog Number	Length (cm)	Outer Diameter (mm)
	w/o spike	XCRF 40622	230	2.4
	fenestrated			
	ellipsoid cups			
	w/spike	XCRF40722	230	2.4
	fenestrated			
	ellipsoid cups			1/box

Bronchoscope Biopsy Forceps

	Description	Catalog No.	Length (cm)	Outer Diameter (mm)
	w/o spike	XCRF 40115	160	1.7
	fenestrated			
	ellipsoid cups			1/box

Cold Biopsy Forceps - Coated Sheath

	Description	Catalog No.	Length (cm)	Outer Diameter (mm)
	w/out spike	XCRF 40622C	230	2.4
	fenestrated			
	ellipsoid cups			
	coated metal spiral			
	w/spike	XCRF 40722C	230	2.4
	fenestrated			
	ellipsoid cups			1/box
	coated metal spiral			

Hot Biopsy Forceps

	Description	Catalog No.	Length (cm)	Outer Diameter (mm)
	w/o spike	XHRF 40622	230	2.4
	fenestrated			
	ellipsoid cups			1/box

Retrieval Devices

»» Reusables & Single Use
»» Steel Wires and Baskets

»» Autoclavable *
»» Durable Handle *

* Reusable Versions

Tripod Graspers



Description	Catalog No.	Length (cm)	Outer Diameter (mm)
Single Use	XFBG-01-23-150	150	2.3
Single Use	XFBG-01-24-220	220	2.3
Single Use	XFBG-02-17-150	150	1.7
10/box			
Reusable	XRFBG-01-23-150	150	2.3
Reusable	XRFBG-01-23-220	220	2.3
Reusable	XRFBG-02-17-150	150	1.7
1/box			

Foreign Body Baskets - Elliptical



Description	Catalog No.	Length (cm)	Outer Diameter (mm)
Single Use	XFBBE 01-23-150	150	2.3
Single Use	XFBBE 01-23-220	220	2.3
Single Use	XFBBE 02-17-150	150	1.7
10/box			
Reusable	XRFBBE 01-23-150	150	2.3
Reusable	XRFBBE 01-23-220	220	2.3
Reusable	XRFBBE 02-17-150	150	1.7
1/box			

Specify: LARGE=30mm or MEDIUM=22mm Basket Size

Foreign Body Baskets - Hexagonal



Description	Catalog No.	Length (cm)	Outer Diameter (mm)
Single Use	XFBBH 01-23-150	150	2.3
Single Use	XFBBH 01-23-220	220	2.3
Single Use	XFBBH 02-17-150	150	1.7
10/box			
Reusable	XRFBBH 01-23-150	150	2.3
Reusable	XRFBBH 01-23-220	220	2.3
Reusable	XRFBBH 02-17-150	150	1.7
1/box			

Specify: LARGE=30mm or MEDIUM=22mm Basket Size

Hemo-Probe

Single-Use Bipolar Coagulation Probe



- »» Round Distal Tip
- »» Spiral Electrode
- »» Hydrophilic Distal Section
- »» Irrigation Port
- »» Non-Stick Ceramic Tip
- »» Compatible w/ Various Generators

Bipolar Coagulation Probe

	Catalog No.	Outer Diameter (FR)	mm	Working Length (cm)	Plug Type
	EL 712.350	7	2.3	350	Single
	EL 1013.350	10	3.3	350	Single
	EL722.350	7	2.3	350	Dual
	EL 1023.350	10	3.3	350	Dual
					1/box

Please specify Generator Type upon placing your order: Olympus, Valleylab®, Everest Medical BiCoag®, Bicap®, Bicap II, Bard®, Microvasive®, Pentax, etc.

Polyp Catch

- »» Two Removable Chambers
- »» Easy to Use
- »» Very Economical



Polyp Catches – Single Use

Description	Catalog No.
Two-Chamber Polyp Catch	EB12555
Removable Chambers	12/bag



Injection/Sclerotherapy Needles

- »» Ergonomic Handle
- »» Lateral Irrigation Port



Supplied Sterile


Single Use Injection Needles

Description	Catalog No.	Needle Gauge/Length	Working Length (cm)	Sheath Diameter (mm)
w/irrigation port	SN-01-22-230	22Ga/5mm	230	2.3

10/box

Transbronchial Aspiration Needles

- »» Atraumatic, Rounded Tip
- »» Various Needle Diameters
- »» Needle Extension Lock



Supplied Sterile

Single Use T.B.A. Needles

Description	Catalog No.	Needle Gauge	Working Length (cm)	Diameter (mm)
TBA Needle	ANSO-19-W	19	140	1.7
TBA Needle	ANSO-20W	20	140	1.7
TBA Needle	ANSO-21-W	21	140	1.7
TBA Needle	ANSO-22-W	22	140	1.7

5/box

Single Use Polypectomy Snares



- »» Three Ring Handle
- »» Braided Wire Loops

- »» Single Use
- »» Olympus Hot Plug

Single Use - 2.4mm, 220cm

Description		Catalog No.	Loop Diameter (mm)	Length (cm)	Sheath Diam. (mm)
	Oval Mini Loop	XPPS-710	10	220	2.4
	Oval Loop	XPPS-715	15	220	2.4
	Oval Loop	XPPS-725	25	220	2.4
	Oval Loop	XPPS-730	30	220	2.4

10/box

Single Use - 1.7mm, 150cm

Description		Catalog No.	Loop Diameter (mm)	Length (cm)	Sheath Diam. (mm)
	Oval Mini Loop	XPPS-810	10	150	1.7
	Oval Loop	XPPS-815	15	150	1.7
	Oval Loop	XPPS-825	25	150	1.7
	Oval Loop	XPPS-830	30	150	1.7

10/box

Single Use - Rotatable

Rotatable Oval Loop					
	Description	Catalog No.	Loop Diameter (mm)	Length (cm)	Sheath Diam. (mm)
	Mini Oval Loop	HIC-SODRT-23	13	220	2.5
	Oval Loop	HIC-SODRT-24	25	220	2.5
	Oval Loop	HIC-SODRT-25	32	220	2.5

10/box

Reusable Polypectomy Snares



- »» Braided Wire Loops
- »» Olympus Hot Plug

- »» Autoclavable
- »» With Flush Port

Reusable - 2.4mm, 220cm

Description	Catalog No.	Loop Diameter (mm)	Length (cm)	Sheath Diam. (mm)
	Mini Oval Loop	XRPPS-71023	10	220
	Oval Loop	XRPPS-71523	15	220
	Oval Loop	XRPPS-72523	25	220
	Oval Loop	XRPPS-73023	30	220

1/box

Reusable - 1.7mm, 150cm

Description	Catalog No.	Loop Diameter (mm)	Length (cm)	Sheath Diam. (mm)
	Mini Oval Loop	XRPPS-71017	10	150
	Oval Loop	XRPPS-71517	15	150
	Oval Loop	XRPPS-72517	25	150
	Oval Loop	XRPPS-73017	30	150


1/box

Channel Cleaning Brushes

- »» Stainless Steel Inner Wire
- »» Metal Spiral or Nylon Sheath

- »» Sturdy Hand Loop
- »» Reusable & Single-Use

Channel Cleaning Brushes - Reusable

	Description	Catalog No.	Sheath	Diameter	Brush Diameter	Working Length
	Metal Spiral	CBS3220		2.5 mm	7 mm	220 cm
	Sheath	CBS2150		2.5 mm	5 mm	150 cm
		CBS2220		2.5 mm	5 mm	220 cm
		CBS1150		1.7 mm	3 mm	150 cm
	Teflon®	CBN3220		2.5 mm	7 mm	220 cm
	Coated Sheath	CBN2150		2.5 mm	5 mm	150 cm
		CBN2220		2.5 mm	5 mm	220 cm
		CBN1150		1.7 mm	3 mm	150 cm
1/box						

Channel Cleaning Brushes - Single-Use

	Description	Catalog No.	Sheath	Diameter	Brush Diameter	Working Length
	Plastic	CBD3230		1.8mm	6.0mm	230 cm
	Sheath					
	Single Head					
	Plastic	CBDD3240		1.8mm	6.0mm	230 cm
	Sheath					
	Dual Head					
10/box						

Cytology Brushes

Single-Use Cytology Brushes - 3mm Diameter

	Description	Catalog No.	Application/ Length cm	Sheath Diameter	Working Length
	3mm Brush	CTB2150	Gastro/ 150cm	2.5 mm	150 cm
	Diameter,	CTB2220	Colon/ 220cm	2.5 mm	220 cm
	Stainless Steel	CTB1110	Bronch/ 110cm	1.7 mm	110 cm
	Wire with	CTB1150	Pedi-Gastro/ 150cm	1.7 mm	150 cm
	Teflon® Outer Sheath				
10/box					

Patient Accessories



Privacy Pants for Colonoscopy Procedures - Single Use

	Description	Size	Catalog No.
	Divided Rear Panel for Easy Access.	One Size	KOL-U00-A1
	Breatheable, water-resistant, triple-layer fabric.	Fits Most	
	Bacteriostatic.		
	Non-sterile.		
	Endo-Bibs	(available outside the U.S. only)	GAS-U00-01

10/box

Bite Blocks - Single Use

	Description	Opening Diameter (mm)	Strap Width (inch)	Maximum Dilator (Fr)	Catalog No.
	Polypropylene	26 X 20	5/8	18	EB61769
	Bite Block				
	with Latex-Free Strap				24/bag